

SPECIFICATIONS

Waveform Specifications



IFC (Interferential 4 pole)

Interferential Current is a medium frequency waveform. Current comes out of two channels (four electrodes). The currents cross each other in the body area that requires treatment. The two currents interfere with each other at this crossing point, resulting in a modulation of the intensity (the current intensity increases and decreases at a regular frequency).

Output Mode:	Pads
Carrier Frequency	Fixed
Beat Frequency	0-200 Hz
Sweep Time	15 seconds
Sweep Low Beat Frequency	1-200 Hz
Sweep High Beat Frequency	1-200 Hz
Scan Percentage	Static, 10%, 40%, 100%
Vector Scan	Manual (0°-90°), Automatic (40% and 100%), Off
Amplitude	0-100 mA RMS into 500 ohm
Maximum Treatment Time	1-60 Minutes
Available on Channel	1&2, 3&4



IFC (Interferential 2 pole) Premodulated

Premodulated Current is a medium frequency waveform. Current comes out of one channel (two electrodes). The current intensity is modulated: it increases and decreases at a regular frequency (the Amplitude Modulation Frequency).

Output Mode:	Pads
Output Intensity	0-100 mA
Carrier Frequency	Fixed
Beat Fixed (Sweep Off)	1-200 Hz
Sweep Low Beat Frequency	1-149 Hz
Sweep High Beat Frequency	81-200 Hz
Cycle Time	Continuous, 5/5, 4/12, 10/10, 10/20, 10/30, 10/50
Mode Selection	CC or CV*
Treatment Time	1-60 Minutes
Available on Channel	1&2, 3&4

*CC= Constant Current
CV= Constant Voltage

SPECIFICATIONS

TENS- Symmetrical Biphasic

The Symmetrical Biphasic waveform has a short pulse duration and is capable of strong stimulation of nerve fibers in the skin and muscle. This waveform is often used in portable muscle stimulation units and some TENS devices. Because of its short pulse duration, the patient typically tolerates the current well, even at relatively high intensities.

Output Mode Pads
Output Intensity 0-80 mA
Phase Duration Adjustable 20-1,000 μ sec
Frequency 1-250 Hz
Mode Selection CC or CV*
Burst Frequency 0-4 bps
Frequency Modulation 0-250 Hz
Amplitude Modulation Off, 40%, 60%, 80% and 100%
Treatment Time 1-60 minutes

DANGER



- Stimulus delivered by the TENS waveforms of this device, in certain configurations, will deliver a charge of 25 microcoulombs (μ C) or greater per pulse and may be sufficient to cause electrocution. Electrical current of this magnitude must not flow through the thorax because it may cause cardiac arrhythmia.



High Voltage Pulsed Current (HVPC)

The High Volt waveform has a very brief pulse duration characterized by 2 distinct peaks delivered at high voltage. The waveform is monophasic (current flows in one direction only). The high voltage causes a decreased skin resistance making the current comfortable and easy to tolerate.

Output Intensity 0-500 V
Output Mode Pads or Probe
Polarity Positive or Negative
Ramp 0.5 sec, 1 sec, 2 sec, 5 sec
Display Peak Current or Volts
Sweep Continuous, 80/120 pps, 1/120 pps, 1/10 pps
Frequency 10-120 pps
Cycle Time	... 5/5, 4/12, 10/10, 10/20, 10/30, 10/50, Continuous
Treatment Time 1-60 Minutes
Available on Channels 1, 2, 3, 4

*CC= Constant Current
CV= Constant Voltage

SPECIFICATIONS

Russian

Russian Current is a sinusoidal waveform, delivered in bursts or series of pulses. This method was claimed by its author (Kots) to produce maximal muscle strengthening effects without significant discomfort to the patient.

Output Intensity	0-100 mA
Output Mode	Pads
Channel Mode	Single, Reciprocal, Co-Contract
Duty Cycle	10%, 20%, 30%, 40%, 50%
Mode Selection	CC or CV*
Anti-Fatigue	Off or On
Cycle Time	5/5, 4/12, 10/10, 10/20, 10/30, 10/50, Continuous
Burst Frequency (Anti-Fatigue Off)	20-100 pps
Ramp	0.5, 1, 2 and 5 seconds
Treatment Time	1-60 minutes
Available on Channels	1, 2, 3, 4

Microcurrent

Microcurrent is a monophasic waveform of very low intensity that closely simulates the electrical current generated by the human body. Microcurrent can be applied via electrodes or probe.

Output Intensity	0-1000.0 μ A
Output Mode	Pads or Probe
Polarity	Positive, Negative or Alternating
Treatment Time	1-60 Minutes
Available on channels	1, 2, 3, 4

*CC= Constant Current
CV= Constant Voltage

SPECIFICATIONS

Ultrasound Specifications

Frequency 1.0 MHz, +/- 5%; 3.3MHz, +/- 5%
Duty Cycles 10%, 20%, 50%, Continuous
Pulse Frequency 100Hz
Pulse Duration 1mSec, +/-20%; 2mSec, +/-20%
5mSec, +/-20%

Output Power

5cm² Crystal 0-10 Watts at 1 & 3.3 MHz
10cm² Crystal 0-20 Watts at 1MHz
0-10 Watts at 3.3 MHz
2cm² Crystal 0-4 Watts at 1 & 3.3 MHz
1cm² Crystal 0-2 Watts 3.3 MHz
Amplitude 0 to 2.5 W/cm² in continuous mode,
0-3 w/cm² in pulsed modes
Output accuracy ± 20% above 10% of maximum
Temporal Peak to Average Ratios:
2:1, ± 20%, at 50% Duty Cycle
5:1, ± 20%, at 20% Duty Cycle
9:1, ± 20%, at 10% Duty Cycle
Beam Nonuniformity Ratio 5.0 : 1 maximum
Beam Type Collimating
Effective Radiating Areas. 10cm² Crystal - 8.5cm², +/- 1.5
5cm² Crystal - 4.0cm², +/- 1.0
2cm² Crystal - 1.8cm², +0.2/-0.4
1cm² Crystal - 0.8cm², +0.2/-0.4
Treatment Time 1-30 Minutes