


SURGICAL SPECIALTIES™

Bone Wax

Softer and easier to handle, form and apply than ever before.

The peel-open foil packet is specifically designed for easy dispensing and provides optimum safety and convenience for use in the operating room.

See for yourself how Surgical Specialties™ Bone Wax compares to Ethicon™ Bone Wax.



Surgical Specialties™ Bone Wax

Ethicon™ Bone Wax

Description

A sterile mixture of beeswax, paraffin, and isopropyl palmitate, a wax-softening agent. It has a waxy odor.



A sterile mixture of beeswax, paraffin, and isopropyl palmitate, a wax-softening agent. It is opaque and has a waxy odor.

Action/Performance

Achieves local hemostasis of bone by acting as a mechanical sealant (tamponade). It does not act biochemically and is minimally resorbable.



Achieves local hemostasis of bone by acting as a mechanical (tamponade) barrier. It does not act biochemically and is minimally resorbable.

Indications

For control of bleeding from bone surfaces during surgical operations.



May be used for the control of bleeding from bone surfaces.

Turn over for more info ➞

Directions & Dosage

Use bone wax immediately after removal from the package. Using aseptic technique, take 1 to 1.5 grams of bone wax, manipulate with the fingers to soften, and apply to bone surface. Optimum working temperature: 70° - 74° F (21° - 23° C).




Bone wax should be used immediately after removal from the package. Using aseptic technique, bone wax should be warmed to desired consistency by manipulation with the fingers or by immersion of the unopened foil packet in a warm sterile solution.

How Supplied

Available sterile in bar form in individual foil packets, 2.5 grams each.



Available sterile in individual foil envelopes, each containing 2.5 grams, and packaged in an individually sealed overwrap packet.



Item #	Description	Packaged
903	White Bone Wax, 2.5 grams each	12 sterile peel-open foil packets per dispenser box
903B	White Bone Wax, 2.5 grams each	Bagged sterile foil packets (your choice of quantity)

To learn more about Surgical Specialties™ Bone Wax or to place an order, contact your local sales representative or our customer service team at:

877.991.1110 (USA and Canada only), **+1.610.404.3520** (International),

or email sscbonewax@surgicalspecialties.com