




Surgical Stainless Steel Suture

Conventional Cutting - Sternum


Needle	Material	Length	none	none	4	5	6	7
CCS 48mm 1/2 Circle 		4-18" ~ 45cm	M652G33	M653G33 M654G33 M655G33				
		2-18" ~ 45cm		M657G33				
	Monofilament B & S 22	4-18" ~ 45cm			M652G			
	Monofilament B & S 20	4-18" ~ 45cm				M653G		
	Monofilament B & S 20	2-18" ~ 45cm				M657G		
	Monofilament B & S 19	4-18" ~ 45cm					M654G	
	Monofilament B & S 18	4-18" ~ 45cm						M655G


Needle	Material	Length	none	5
CCS-1 60mm 1/2 Circle 		2-30" ~ 75cm	M436G33	
	Monofilament B & S 20	2-30" ~ 75cm		M436G


Reverse Cutting


Needle	Material	Length	none	3-0
CP-1 36mm 1/2 Circle 		20" ~ 50cm	497G33	
	Monofilament B & S 30	20" ~ 50cm		497G

Sutures : Non Absorbable / Surgical Stainless Steel Suture

Needle	Material	Length	3-0
FS-1 24mm 3/8 Circle 	Monofilament B & S 30	18" ~ 45cm	631G


Needle	Material	Length	none	3-0	4-0
FS-2 19mm 3/8 Circle 		18" ~ 45cm	603G33		
	Monofilament B & S 30	18" ~ 45cm		607G	
	Monofilament B & S 32	18" ~ 45cm			603G


Needle	Material	Length	none	3-0
FSL 30mm 3/8 Circle 		18" ~ 45cm	606G33	
	Monofilament B & S 30	18" ~ 45cm		606G

Needle	Material	Length	none	2-0
FSLX 40mm 3/8 Circle 		18" ~ 45cm	618G33	
	Monofilament B & S 28	18" ~ 45cm		618G


Needle	Material	Length	2

Sutures : Non Absorbable / Surgical Stainless Steel Suture


Needle	Material	Length	2	
LR 75mm LR 3/8 Circle 	Monofilament B & S 24	30" ~ 75cm	467G	

Needle	Material	Length	none	2
LS-1 45mm 3/8 Circle 		4-20" ~ 50cm	M400G33	
	Monofilament B & S 24	4-20" ~ 50cm		M400G


Straight Cutting Needles


Needle	Material	Length	5-0	
KS 60mm Straight 	Monofilament B & S 35	8-16" ~ 40cm	M645G	

SUTUPAK Pre-Cut Sutures


Needle	Material	Length	0	2-0	3-0
SUTUPAK Pre-Cut Sutures N/A 	Monofilament B & S 26	17-18" ~ 45cm	SS26G		
	Monofilament B & S 28	17-18" ~ 45cm		SS28G	
	Monofilament B & S 30	17-18" ~ 45cm			SS30G


Taper Point

Needle	Material	Length	2-0
CT-1 36mm 1/2 Circle 	Monofilament B & S 28	8-18" ~ 45cm	M682G


Needle	Material	Length	2-0
CTX 48mm 1/2 Circle 	Monofilament B & S 28	8-24" ~ 60cm	M414G


TAPERCUT - Straight Cutting


Needle	Material	Length	4-0
TS_TS 48mm TS Straight  <p>Polypropylene</p>	Monofilament B & S 32	2-14" ~ 35cm	540G ⁴¹

Needle	Material	Length	4-0
V-26 26mm 3/8 Circle  <p>Polypropylene</p>	Monofilament B & S 32	2-14" ~ 35cm	539G ⁴¹

TAPERCUT

Needle	Material	Length	1
V-37 40mm 1/2 Circle 	Monofilament B & S 25	4-18" ~ 45cm	M660G


Needle	Material	Length	4	5	6
V-40 48mm 1/2 Circle 	Monofilament B & S 22	4-18" ~ 45cm	M651G		
	Monofilament B & S 20	4-18" ~ 45cm		M650G	
	Monofilament B & S 20	2-18" ~ 45cm		M656G	
	Monofilament B & S 19	4-18" ~ 45cm			M649G


Needle	Material	Length	5
V-56 60mm 1/2 Circle 	Monofilament B & S 20	4-30" ~ 75cm	M435G


Temporary Pacing Wire

Needle	Material	Length	2-0
--------	----------	--------	-----

Sutures : Non Absorbable / Surgical Stainless Steel Suture


Needle	Material	Length	2-0
BB 17mm SKS 60mm 3/8 Circle 	Dk. Blue Braided Stainless Steel	24" ~ 60cm	TPW20 ¹


Needle	Material	Length	2-0
RB-1 17mm SKS 60mm 1/2 Circle Straight Cutting 	Dk. Blue Braided Stainless Steel	24" ~ 60cm	TPW10 ¹

Needle	Material	Length	0	2-0
SH 26mm SKS-3 90mm 1/2 Circle Straight Cutting 	Dk. Blue Braided Stainless Steel	24" ~ 60cm	TPW50 ¹	TPW90
	Lt. Blue & Dk. Blue Braided Stainless Steel	2-24" ~ 60cm	TPW52	TPW92
	Dk. Blue Braided Stainless Steel	24" ~ 60cm		TPW30 TPW70
	Lt. Blue & Dk. Blue Braided Stainless Steel	2-24" ~ 60cm		TPW32 TPW72

Needle	Material	Length	2-0
--------	----------	--------	-----

Sutures : Non Absorbable / Surgical Stainless Steel Suture

Needle	Material	Length	2-0
BB-1 22mm SKS-3 90mm 3/8 Circle 	Dk. Blue Braided Stainless Steel	24" ~ 60cm	TPW40 ^{1 68}
	Lt. Blue & Dk. Blue Braided Stainless Steel	2-24" ~ 60cm	TPW42

Needle	Material	Length	0	2-0
TE 31mm SKS-3 90mm 3/8 Circle 	Lt. Blue & Dk. Blue Braided Stainless Steel	2-24" ~ 60cm	TPW62 ^{1 11}	
	Dk. Blue Braided Stainless Steel	24" ~ 60cm		TPW80

¹ One dozen packets per box

¹¹ Multistrand 2" ~ 60 cm

⁴¹ Pull out wire set (includes one polypropylene button)

⁶⁸ 1.5" long scoremarks; others .875"