

# BD Disposal Solutions

## Product Catalog



Helping all people  
live healthy lives

# The BD Recykleen™ Story

Being sensitive to environmental issues is very important for all people and businesses. BD recognizes this need and is searching for ways to be a better environmental steward while helping our customers. One result is the BD Recykleen product line.

These products are made from post-industrial or post-consumer recycled plastic. The amount of recycled material depends upon the product. While BD strives to use as much recycled plastic as possible in each product, there are practical things to consider like:

- **Supply** – Availability of recycled resin sufficient to meet the market demand
- **Quality** – Consistency and reliability of recycled resin material
- **Performance** – Products made with recycled resin must meet the same stringent performance requirements as those from virgin plastic
- **Visual** – Some products need transparency, which cannot be provided with most recycled plastics

This product line will continue to grow and includes:



Collector not included.

**BD Recykleen™  
Cabinets**



**BD RCRA  
Collectors**



**BD Recykleen™  
L and XL Collectors**



**BD Recykleen™  
Trolleys**







**BD Recykleen™  
Patient Room**

# TABLE OF CONTENTS

Waste Stream Segregation .....	4
How to Choose the Right Collector .....	5
Patient Room .....	6
Cabinets and Accessories .....	7
Large and Extra Large Sharps Collectors .....	8
Trolleys .....	9
Multi-use Nestable .....	10
Multi-use One Piece .....	11
Phlebotomy .....	11
Pharmaceutical .....	12
Chemotherapy .....	13
RCRA DOT Collectors .....	14
RCRA Non DOT Collectors .....	15

Look for these symbols throughout the catalog to easily identify key features and benefits:




# WASTE STREAM SEGREGATION

BD Disposal Solutions offers collectors to meet the disposal needs of the four categories of waste generated in medical settings. There are two main categories of waste: regulated medical waste and hazardous waste. Regulated medical waste includes Sharps / Infectious Waste, Pharmaceutical (non-hazardous) and Chemotherapy waste.

The hazardous waste category is for the collection of RCRA (Resource Conservation and Recovery Act by the U.S. EPA) waste.

This table indicates general guidelines for proper waste stream segregation.



# HOW TO CHOOSE THE RIGHT COLLECTOR

Waste Stream	Type	Capacity	Available in BD Recykleen™	Access Opening Size (inches)	Sample Products that fit in Collector
Sharps and Infectious Waste	Patient Room	5.4 qt with Counterbalanced Door	●	2.0 X 8.0	Blood collection tubes and collection assembly, syringes up to 60 mL with needle attached
		2 gal with Counterbalanced Door		2.0 X 8.0	
		3 gal with Counterbalanced Door	●	2.0 X 8.0	
		5.4 qt side entry		1.5 X 5.2	Syringes and slides
	Phlebotomy	1 qt		1.2 X 2.5	1 mL - 10 mL syringes, blood collection tubes and assembly, slides, pipettes
		1.5 qt		1.6 X 1.6	
	Multi-use Nestable	8 qt		1.9 Dia Open (305417 only); 2.3 Dia Petals	1 mL - 20 mL syringes with needle, blood collection tubes, slides, pipettes
		14 qt		2.6 Dia Open (305465); 4.4 Dia Petals (305464); 4.5 Dia Open (305456)	1 mL - 60 mL syringes with needle, blood collection assembly, blood collection tubes, slides, pipettes
		6 gal		4.5 Dia Open (305457); Petals (305465)	All syringe sizes with needle attached, slides, urine collection assembly, blood culture bottles, blood collection assembly
	Multi-use One Piece	1.4 qt		1.8 Dia Petals	Small syringes, needles, pipettes, slides
		3.3 qt		1.8 Dia Petals	
		6.9 qt		1.8 Dia Petals	
		8.2 qt		1.8 Dia Petals	
	Large and Extra Large	5 gal		4.2 Dia Open	All syringe sizes with needle attached, slides, urine collection assembly, blood culture bottles, blood collection assembly
		9 gal Hinge Top	●	9.0 X 10.0	All items mentioned above plus small laparoscopic instruments, phlebotomy procedure bottles, dialysis catheters
		9 gal Slide Top	●	5.0 X 7.5	
		17 gal Hinge Top	●	10.5 X 12.5	Medical devices, spinal tap devices, laparoscopy instrumentation
		19 gal Slide Top	●	6.0 X 8.5	
Pharmaceutical	Patient Room	5.4 qt with Counterbalanced Door	●	2.0 X 8.0	All syringe sizes with needle attached, small IV bags 50 mL - 250 mL, medication vials, ampules
		3 Gal with Counterbalanced Door	●	2.0 X 8.0	All syringe sizes with needle attached, small IV bags 50 mL - 250 mL
	Multi-use One Piece	3 gal		4.2 Dia Open	All syringe sizes and bags up to 250 mL
	Large and Extra Large	9 gal Hinge Top		9.0 X 10.0	All syringe sizes, all IV bag sizes with tubing attached, all medication vials, ampules, TPN bags, IV glass bottles
		9 gal Slide Top		5.0 X 7.5	
		17 gal Hinge Top		10.5 X 12.5	
		19 gal Slide Top	●	6.0 X 8.5	
Chemotherapy	Multi-use One Piece	3 gal		4.2 Dia Open	All syringe sizes with needle attached, slides, urine collection assembly, blood culture bottles, blood collection assembly
		5 gal		4.2 Dia Open	
		5 gal Screw Cap		3.8 Dia Open	
	Large and Extra Large	9 gal Hinge Top		9.0 X 10.0	All syringe sizes, all IV bag sizes with tubing attached, all medication vials, ampules, TPN bags, IV glass bottles
		9 gal Slide Top		5.0 X 7.5	
		17 gal Hinge Top		10.5 X 12.5	
RCRA	Multi-use One Piece	3 gal		4.2 Dia Open	Smaller size of solid RCRA hazardous materials plus liquids if contained in a leak resistant container
		5 gal		4.2 Dia Open	
	Large and Extra Large	9 gal Hinge Top		9.0 X 10.0	Solid RCRA hazardous materials plus liquids if contained in a leak resistant container
		9 gal Slide Top		5.0 X 7.5	
		17 gal Hinge Top		10.5 X 12.5	
		19 gal Slide Top		6.0 X 8.5	
	PGII	13 gal		6.0 Dia Open	Solid RCRA hazardous materials plus liquids if contained in a leak resistant container
		16 gal		6.0 Dia Open	



# PATIENT ROOM

BD Patient Room Collectors are available in BD Recykleen™ and BD standard brands and come in a variety of sizes and colors. All BD Patient Room Collectors have visible fill level indicators and handles for safe transportation.

## Key product features include:

- Compatible with locking cabinets to provide security and prevent unwanted access to collector contents
- Counterbalanced door facilitates hands-free disposal and is designed to reduce the potential of sharps injuries
- Designed to improve tamper resistance by preventing hand access to the contents
- Available in BD Recykleen brand



305517



305053



## Patient / Exam Room

Reference #	Product Size	Product Dimensions (H x W x D inches)	Qty Case	Compatible with
<b>BD Recykleen™ Sharps Collectors</b>				
305056	5.4 qt BD Recykleen Pearl Counterbalanced Door	10.75 x 12 x 4.5	20	305096, 305098, 305017, 305447, 305021
305057	5.4 qt BD Recykleen Red Counterbalanced Door	10.75 x 12 x 4.5	20	
305053	3 gal BD Recykleen Pearl Counterbalanced Door	15.75 x 13.5 x 6	10	305097, 305099, 305420 305409, 305588, 305552
305054	3 gal BD Recykleen Red Counterbalanced Door	15.75 x 13.5 x 6	10	
305517	5.4 qt Red Counterbalanced Door	10.75 x 12 x 4.5	20	305096, 305098, 305017, 305447, 305021
305551	5.4 qt Clear Counterbalanced Door	10.75 x 12 x 4.5	20	
305543	5.4 qt Pearl Counterbalanced Door	10.75 x 12 x 4.5	20	
305426	5.4 qt Red Counterbalanced Door	10.75 x 12 x 4.5	12	
305427	5.4 qt Clear Counterbalanced Door	10.75 x 12 x 4.5	12	
305428	5.4 qt Pearl Counterbalanced Door	10.75 x 12 x 4.5	12	
305422	2 gal Pearl Counterbalanced Door	11.5 x 13.5 x 6	10	305097, 305099, 305420, 305409, 305588, 305552
305435	2 gal Red Counterbalanced Door	11.5 x 13.5 x 6	10	
305423	3 gal Pearl Counterbalanced Door	15.75 x 13.5 x 6	10	
305436	3 gal Red Counterbalanced Door	15.75 x 13.5 x 6	10	
305443	5.4 qt Red Side Entry	10.75 x 12 x 4.5	20	305447, 305021
305444	5.4 qt Pearl Side Entry	10.75 x 12 x 4.5	20	
305425	5.4 qt Pearl Side Entry	10.75 x 12 x 4.5	12	

See Accessories on page 7

# CABINETS AND ACCESSORIES

BD Recykleen™ Cabinets are available in two sizes: small for use with 5.4 quart collectors and large for use with 2 and 3 gallon collectors. Each size is available in both green and brown semi-transparent colors.



## Key product features include:

- A universal lock providing security and preventing unwanted access to collectors
- A one-piece transparent door facilitating easy cleaning
- The ability to clearly see the collector fill level indication
- A two-year warranty providing product replacement in case of manufacturer defects



## BD Recykleen™ Cabinets

Reference #	Product Size	Product Dimensions ( H x W x D inches )	Qty Case	Compatible with
305096	5.4 qt BD Recykleen Brown Cabinet	14 x 6 x 15	1	305056, 305057, 305517, 305551, 305543, 305426, 305427, 305428
305098	5.4 qt BD Recykleen Green Cabinet	14 x 6 x 15	1	
305097	2 and 3 gal BD Recykleen Brown Cabinet	15 x 7 x 19	1	305053, 305054, 305422, 305435, 305423, 305436
305099	2 and 3 gal BD Recykleen Green Cabinet	15 x 7 x 19	1	

## Brackets

Brackets are designed to secure the collector to a wall or cart. Some brackets include additional security via a locking mechanism.

## Stabilizers

Stabilizers are designed to stabilize collectors on a countertop or other flat surface.

## Brackets and Stabilizers


Reference #	Product Size	Product Dimensions ( H x W x D inches )	Qty Case	Compatible with
305447	Locking Bracket	9.5 x 9.5 x 4.25	10	All 5.4 qt
305409	Locking Bracket	4.25 x 14 x 4.75	10	2 and 3 gal Counterbalanced Door
305588	Wire Bracket	6 x 12.5 x 6	6	2 and 3 gal Counterbalanced Door
305021	Stabilizer	3 x 13 x 6.5	5	All 5.4 qt
305552	Stabilizer	3 x 13.5 x 7.5	5	2 and 3 gal Counterbalanced Door
305448	Glove Box	7 x 12 x 4.25	6	Attaches to Cabinets and 5.4 qt Locking Bracket
305495	XL Security Bracket	2 x 10.5 x 8	12	Multi-use One Piece 5 gal, 3 and 5 gal (Plug Cap)
305497	Quick Release Strap	1.5 x 8 x 1	12	Multi-use One Piece 3.3 qt, 6.9 qt, 8.2 qt, 3 gal (Plug Cap), 5 gal

# LARGE AND EXTRA LARGE SHARPS COLLECTORS

BD Large and Extra Large Collectors are available in BD Recykleen™ and BD standard brands. These collectors are fully compatible with trolleys to provide stability and mobility (see page 9 for details).



305023



## Key product features include:

- Temporary and final closure features
- Visual, audible, tactile (VAT) perimeter locks to ensure that the lid is securely attached to the base and to provide security against tampering
- Double wall handles for safe transportation
- Clear tops to ensure visibility of collector fill level
- Easy accessibility for disposal of large and bulky waste
- Dual Opening Hinge Top designed so that the small or the large side can remain open and the other locked based on the size of access opening desired



305616

## Large and Extra Large Sharps

Reference #	Product Size	Product Dimensions ( H x W x D inches )	Qty Case	Compatible with
<b>BD Recykleen™ Sharps Collectors</b>				
305024	9 gal BD Recykleen Hinge Top	18.75 x 17.75 x 11.75	8	305090
305025	9 gal BD Recykleen Slide Top	18.75 x 17.75 x 11.75	8	305090, 305091
305023	17 gal BD Recykleen Hinge Top	23.5 x 20 x 14.75	5	305092
305026	19 gal BD Recykleen Slide Top	26.25 x 20 x 14.75	5	305092, 305093
305615	9 gal Clear Hinge Top	18.75 x 17.75 x 11.75	8	305090
305601*	9 gal Clear Hinge Top with gasket	18.75 x 17.75 x 11.75	8	
305616	9 gal Clear Slide Top	18.75 x 17.75 x 11.75	8	305090, 305091
305602*	9 gal Clear Slide Top with gasket	18.75 x 17.75 x 11.75	8	
305610	17 gal Clear Hinge Top	23.5 x 20 x 14.75	5	305092
305665*	17 gal Clear Hinge Top with gasket	23.5 x 20 x 14.75	5	
305609	19 gal Clear Slide Top	26.25 x 20 x 14.75	5	305092, 305093
305666*	19 gal Clear Slide Top with gasket	26.25 x 20 x 14.75	5	

\* Meets or exceeds **US DOT** packaging requirements 49CFR173.24, 49CFR173.24a, and 49CFR173.134.

See Trolleys on page 9



# TROLLEYS

Our full line of BD Recykleen™ Trolleys are available in basic and foot-operated versions for use with large and extra large collectors. Trolleys provide stability and mobility for large and extra large collectors.



## Key product features include:

- A design that reduces the risk of collectors tipping over
- Wheels that enable the collector to be conveniently placed where it will be most used
- Locking wheels, ensuring the product is mobile only when you want it to be
- The ability to open and close the collector for hands-free disposal (foot-operated trolleys only)
- A two-year warranty providing product replacement in case of manufacturer defects



## BD Recykleen™ Trolleys

Reference #	Product Size	Product Dimensions (H x W x D inches)	Qty Case	Compatible with
305090	9 gal Basic BD Recykleen Trolley	23.5 x 16.3 x 28.6	1	All 9 gal
305091	9 gal Foot-operated BD Recykleen Trolley	23.5 x 16.3 x 28.6	1	All 9 gal slide top
305092	17 and 19 gal Basic BD Recykleen Trolley	24.1 x 18 x 27.5	1	All 17 and 19 gal
305093	19 gal Foot-operated BD Recykleen Trolley	24.3 x 17.7 x 35.2	1	All 19 gal slide top

# MULTI-USE NESTABLE

BD Multi-use Nestable Collectors are available in three sizes and can be used with accessories for added security and stability.

## Key product features include:

- Clear tops with multiple opening styles to accommodate different sizes of waste and environments
- Handles designed for safe transport
- Funnel-style lids designed to prevent accessibility to the collector contents and minimize spillage if the collector tips over
- Open-top lids feature larger diameter openings to accommodate disposal needs in non-patient/ public areas



305457



Regular Funnel  
Clear Top



Large Funnel  
Clear Top



2" diameter Open  
Clear Top



2.5" diameter  
Open Clear Top



4.25" diameter  
Open Clear Top



305417 with  
305570 Stabilizer



305646 with  
305485 Bracket



## Multi-use Nestable

Reference #	Product Size	Product Dimensions (H x W x D inches)	Qty Case	Compatible with
305417	8 qt Open Clear Top	11.25 x 10 x 6.75	24	305485, 305018, 305570
305460	8 qt Regular Funnel Clear Top	11.25 x 10 x 6.75	24	
305463	8 qt Pearl Regular Funnel Clear Top	11.25 x 10 x 6.75	24	
305646	14 qt Small Open Clear Top	11.5 x 12.5 x 8.5	20	
305456	14 qt Open Clear Top	11.5 x 12.5 x 8.5	20	
305464	14 qt Regular Funnel Clear Top	11.5 x 12.5 x 8.5	20	
305480	14 qt Large Funnel Clear Top	11.5 x 12.5 x 8.5	20	
305457	6 gal Open Clear Top	17.5 x 12.5 x 8.5	12	
305465	6 gal Large Funnel Clear Top	17.5 x 12.5 x 8.5	12	

## Brackets and Stabilizers

305485	Bracket	6 x 3.5 x 1.25	12	8 qt, 14 qt and 6 gal
305018	Wire Bracket	9.5 x 10.25 x 7.75	6	All 8 qt
305570	Stabilizer	2 x 10 x 7	5	All 8 qt
305494	Padlock	1.75 x 1 x .25	12	305485, 305018, 305570

See Accessories on page 7

# MULTI-USE ONE PIECE

BD Multi-use One Piece Collectors are available in 5 sizes and can be used with accessories for added security.



305477

## Key product features include:

- Handles designed for safe movement and transportation
- Hinge caps for 6.9 and 8.2 quart collectors designed to provide temporary and final closure features



305557

305488

## Multi-use One Piece

Reference #	Product Size	Product Dimensions (H x W x D inches)	Qty Case	Compatible with
305557	1.4 qt Regular Funnel	7.75 x 3.75 x 3.75	36	
305488	3.3 qt Regular Funnel	5.75 x 8.5 x 5	24	305497
305489	6.9 qt Regular Funnel	11.5 x 8.75 x 5.5	12	
300467*	6.9 qt Regular Funnel	11.5 x 8.75 x 5.5	12	305496
305490	8.2 qt Regular Funnel	13.5 x 8.75 x 5.5	12	305497
300470*	8.2 qt Regular Funnel	13.5 x 8.75 x 5.5	12	
305491	5 gal Open Top	18 x 10.5 x 7.5	8	
300473*	5 gal Open Top	18 x 10.5 x 7.5	8	305495
305477	5 gal Large Funnel	18 x 10.5 x 7.5	8	305497
305577	5 gal Screw Cap	18 x 10.5 x 7.5	8	

## Brackets and Stabilizers

305497	Quick Release Strap	1.5 x 8 x 1	12	3.3 qt, 6.9 qt, 8.2 qt, 5 gal
305496	Security Bracket	1.5 x 8 x 8	12	6.9 qt, 8.2 qt
305495	XL Security Bracket	2 x 10.5 x 8	12	5 gal
305494	Padlock	1.75 x 1 x .25	12	305495, 305496

\* Cannot be autoclaved and not vented

# PHLEBOTOMY

BD Phlebotomy Collectors are available in two sizes and can be used with brackets for added stability.

## Key product features include:

- Built-in needle removal ports designed to prevent contact with sharps during disposal
- Clear tops featuring a clear view of the fill level
- Temporary and final closures for safety and security



305635



305487 with  
305970 Bracket

## Phlebotomy



Reference #	Product Size	Product Dimensions (H x W x D inches)	Qty Case	Compatible with
305635	1 qt Small Open Top	6.5 x 4.75 x 4	60	
305487	1.5 qt Removal Port Top	9 x 4.75 x 4	36	305970
<b>Brackets</b>				
305970	Cup Bracket	3 x 3.5 x 3	12	305635, 305487

# PHARMACEUTICAL

BD Pharmaceutical Collectors feature a variety of lid designs and sizes to collect and dispose of unused, expired or other pharmaceutical waste (see Waste Stream Segregation on pg 4).

## Key product features include:

- Gaskets or liquid absorbing materials to prevent leakage
- Visual, audible, tactile (VAT) perimeter locks to ensure that the lid is securely attached to the base and to add security against tampering
- Handles designed for safe transportation
- Dual Opening Hinge Top designed so that the small or the large side can remain open and the other locked based on the size of access opening desired
- Large and extra large collectors are fully compatible with BD Trolleys (see page 9)
- 5.4 quart and 3 gallon Pharmaceutical Collectors with counterbalanced door:
  - Facilitate hands-free disposal and are designed to reduce the potential for sharps injuries
  - Designed to improve tamper resistance by preventing hand access to the contents
  - Compatible with locking cabinets to provide security and prevent unwanted access to collector contents
  - Available in BD Recykleen brand



## Pharmaceutical

Reference #	Product Size	Product Dimensions ( H x W x D inches )	Qty Case	Compatible with
<b>BD Recykleen™ Sharps Collectors</b>				
305058	5.4 qt BD Recykleen Counterbalanced Door	10.75 x 12 x 4.5	20	305096, 305098, 305017
305055	3 gal BD Recykleen Counterbalanced Door	15.75 x 13.5 x 6	10	305097, 305099, 305420
305330	5.4 qt Counterbalanced Door	10.75 x 12 x 4.5	20	305096, 305098, 305017
305622	3 gal Counterbalanced Door	15.75 x 13.5 x 6	10	305097, 305099, 305420
305633	3 gal Plug Cap	12 x 7.5 x 10.5	12	305495, 305497
305634	9 gal Hinge Top	18.75 x 17.75 x 11.75	8	305090
305329	9 gal Slide Top	18.75 x 17.75 x 11.75	8	305090, 305091
305028	17 gal Hinge Top	23.5 x 20 x 14.75	5	305092
305328	19 gal Slide Top	26.25 x 20 x 14.75	5	305092, 305093

See Accessories on page 7 and Trolleys on page 9

# CHEMOTHERAPY

BD Chemotherapy Collectors feature a variety of lid designs and sizes for disposal of trace chemotherapy waste (see Waste Stream Segregation on page 4). Trace waste is defined by CDC/NIOSH as “those that contain less than 3% by weight of the original quantity of hazardous drugs.”\*

\* CDC/NIOSH Alert Publication No. 2004-165

## Key product features include:

- Gaskets or liquid absorbing materials to prevent leakage
- Visual, audible, tactile (VAT) perimeter locks to ensure that the lid is securely attached to the base and to add security against tampering
- Temporary closure feature and handles designed for safe transportation
- The large and extra large collectors are fully compatible with BD trolleys (see page 9)
- Clear lids with multiple opening styles to accommodate different types of waste and environments
  - Dual Opening Hinge Top designed so that one side can remain open and the other locked based on the size access opening desired





305613



305493



305492



## Chemotherapy

Reference #	Product Size	Product Dimensions (H x W x D inches)	Qty Case	Compatible with
305492	3 gal Plug Cap	12 x 7.5 x 10.5	12	
305493	5 gal Plug Cap	18 x 10.5 x 7.5	8	305495 305497
305578	5 gal One Piece Screw Cap	18 x 10.5 x 7.5	8	
305603	9 gal Clear Hinge Top	18.75 x 17.75 x 11.75	8	305090
305604	9 gal Clear Slide Top	18.75 x 17.75 x 11.75	8	305090, 305091
305614	17 gal Clear Hinge Top	23.5 x 20 x 14.75	5	305092
305613	19 gal Clear Slide Top	26.25 x 20 x 14.75	5	305092, 305093

See Accessories on page 7 and Trolleys on page 9

# RCRA HAZARDOUS WASTE COLLECTORS — DOT


BD RCRA Collectors are intended for Resource Conservation and Recovery Act (RCRA) solid hazardous waste and trace liquids associated with these wastes (see Waste Stream Segregation on page 4).

BD offers two versions of RCRA collectors, **DOT PGII/UN rated** and **non-DOT**, to meet a variety of needs. DOT/UN rated RCRA collectors are designed to be transported without being overpacked by the waste haulers. As most waste haulers will overpack smaller collectors whether they are DOT rated or not, BD offers these in 13 and 16 gallon sizes.

## What is UN?

The UN or *United Nations Recommendations for the Transport of Dangerous Goods Model Regulations* are suggested regulations for controlling the transport of hazardous materials, including medical waste. The regulations define Packing Groups and test requirements to qualify for a UN-compliant package marking.

A UN-compliant package must have a marking that describes its maximum content weight and the Packing Group that it passed.

It looks like:  4H2 / Y26.9 / S / 09  
USA / +CC5871

The regulations describe packaging, labeling and documentation required for transporting hazardous materials, contain a table with codes for each specific type of hazardous material and have several codes for pharmaceutical waste.

## What is DOT?

The DOT or *U.S. Department of Transportation Regulations*, is a government agency regulating fast, safe and efficient transportation of dangerous goods.

## What is PGII?

PG (Packing Group) is a Department of Transportation term that groups hazardous materials according to the danger presented by the materials. There are three packing group classifications.

I- great      II- medium      III- minor

## Key product features of RCRA Hazardous Waste Collectors include:

- Black base with black tops to distinguish the waste stream
- Leak resistant gasketed tops with 7" port opening



## RCRA Hazardous Waste Collectors DOT PGII / UN

Reference #	Product Size	Product Dimensions ( H x W x D inches )	Qty Case
305073	13 gal DOT Port Top	16 x 11 x 21.5	6
305074	16 gal DOT Port Top	16 x 11 x 24.75	6



# RCRA HAZARDOUS WASTE COLLECTORS — Non-DOT

Non-DOT collectors feature a variety of lid designs and sizes and come with a liquid-absorbent material inside the collector for trace waste.\*

\* CDC/NIOSH Alert Publication No. 2004-165

## Key product features of non-DOT Collectors include:



- Made from a minimum of 70% recycled plastic
- Clear tops with multiple opening styles to accommodate different types of waste and environments
  - Dual Opening Hinge Top designed so that one side can remain open and the other locked based on department needs
  - Dual Lock Cap (3 and 5 gallon) designed to provide temporary closure for safe intermediate transportation and secure final closure with two locks
- Visual, audible, tactile (VAT) perimeter locks to ensure that the lid is securely attached to the base and to add security against tampering
- Handles designed for safe transportation
- Temporary closure feature and handles designed for safe transportation
- Large and extra large collectors are fully compatible with BD trolleys (see page 9)



305069



305066



All RCRA collectors must be incinerated for final disposal.

## RCRA Hazardous Waste Collectors

Reference #	Product Size	Product Dimensions (H x W x D inches)	Qty Case	Compatible with
305066	3 gal Plug Cap	12 x 7.5 x 10.5	12	305495
305067	5 gal Plug Cap	18 x 10.5 x 7.5	8	305497
305068	9 gal Clear Hinge Top	18.5 x 17.75 x 11.75	8	305090
305069	9 gal Clear Slide Top	18.5 x 17.75 x 11.75	8	305090, 305091
305071	17 gal Clear Hinge Top	23.5 x 19.5 x 14.25	5	305092
305072	19 gal Clear Slide Top	26.25 x 19.5 x 14.25	5	305092, 305093

See Accessories on page 7 and Trolleys on page 9

## BD ENVIRONMENTAL STEWARDSHIP

BD's approach to corporate citizenship is embedded in our culture, and is intrinsic to our purpose, *"Helping all people live healthy lives."* Responsible corporate citizenship is integrally tied to our global business strategy. Sustained performance and a dedication to improvement are at BD's core. Recognizing that human health and a healthy environment are inseparable, BD is committed to responsible environmental stewardship in our business conduct.




## BD DISPOSAL SOLUTIONS: A COMMITMENT TO QUALITY

BD Medical Waste Collectors are manufactured using polypropylene or polyethylene resins and colorants that are nontoxic and comply with the Toxics in Packaging Clearing House (TPCH) heavy metal limits. No latex materials are used in these products. The collectors are safe for incineration. Under combustion with adequate access of air, the resin burns to water vapor, carbon dioxide and carbon ash.

BD Sharps Collectors meet or exceed the puncture force requirements defined by the American Society for Testing and Materials (ASTM) and British Standards. Sharps collectors are classified as medical devices and are manufactured according to the FDA Quality System Regulation and the International Quality Systems Standard ISO 9001.

BD Sharps Collectors are listed with the FDA and meet pre-market registration requirements in accordance with 21CFR807. BD Sharps Collectors also meet the Occupational Safety and Health Administration (OSHA) Bloodborne Pathogens Standards 29 CFR 1910.1030.



Contact Customer Service at **1-888-237-2762**  
or your BD Sales Representative for additional information



**BD Medical**  
4665 North Avenue  
Oceanside, CA 92056  
[www.bd.com/sharps](http://www.bd.com/sharps)