

Colpo-Master™ Colposcopes

Bovie introduces the new Colpo-Master™ line, a perfect addition to the already popular line of gynecologic products. The line of colposcopes offers superior, 3D, stereoscopic viewing with excellent clarity during gynecologic procedures where millimeters matter. The scopes are available with consistent and dependable zoom range including simple 3 point magnification (3.75x, 7.5x and 15x) or continuous zoom 4 – 27x magnification. For the EMR ready facility the Colpo-Master line is available with Trinocular camera-ready heads to provide quick and easy visual records of treatment. All Colpo-Master Colposcopes are backed by a 7 year limited warranty. Paired with the Aaron 1250-G Electrosurgery System with Smoke Evacuation for a complete GYN LEEP system.

NEW Colpo-Master™ I



Low center of gravity increases stability, while EZ-Roll locking casters provide smooth gliding movement between exam rooms.

CS-103LED includes 3-legged base as shown. Also available as CS-105LED with 5-legged base.

Colpo-Master™ is a registered trademark of United Products and Instruments, Inc; DBA UNICO

COLPO-MASTER™ I ZOOM SUSPENSION ARM COLPOSCOPES

The Colpo-Master™ is perfect for the GYN performing both exams and procedures. The swing-arm design allows easy access to the vaginal canal during surgical procedure.

Colpo-Master™ I features smooth, continuous zoom 4 – 27x magnification with 20x eyepieces and sharp focus for 3-D imaging. Available with Binocular or camera-ready Trinocular heads with 298mm working distance. The field of view (FOV) ranges from 58.5mm at 4x to 8.5mm at 27x. All models with your choice of 3 or 5 leg sturdy mobile bases.

Designed with your viewing comfort in mind; each model has dual diopter adjustments to independently adjust the focus of both Left and Right eyepieces to eliminate eyestrain. Camera-ready Trinocular models allow the user to view through both eyepieces while using the camera.

Colpo-Master™ I's Suspension Arm is rock-steady. To make head positioning easy it comes with locking casters and stainless steel handlebars, standard.

Complete with dust cover, power cord, manual, tool kit and 7 Year Limited Warranty, 1 year on electrical system (excludes bulbs and fuses). Optional accessories are available - please contact your Bovie Representative for additional information.

Item No.	Description
CS-103LED	Colpo-Master™ I Binocular Zoom Suspension-Arm Colposcope
CS-103LED-E	Colpo-Master™ I Binocular Zoom Suspension-Arm Colposcope, 220V
CS-103TLED	Colpo-Master™ I Camera-Ready Zoom Swing-Arm Colposcope
CS-103TLED-E	Colpo-Master™ I Camera-Ready Zoom Swing-Arm Colposcope, 220V
CS-105LED	Same as CS-103LED above, with 5-legged base
CS-105TLED	Same as CS-103TLED above, with 5-legged base

NEW Colpo-Master™ II

C205T-LED Camera-Ready Trinocular head with camera



COLPO-MASTER™ II CENTER-POST MOUNT COLPOSCOPES

Colpo-Master™ II is a compact Center-Post Colposcope with unique, continuous zoom magnification and sharp focusing for crisp 3-D imaging.

Features:

- **HEAD:** Full continuous zoom magnification range: 4 – 27x with 20x eyepieces. Binocular or video-ready Trinocular heads, with 298mm working distance. The field of view (FOV) ranges from 58.5mm at 4x to 8.5mm at 27x
- 3-D Stereoscopic Viewing
- Dual diopters allow user to independently focus each eyepiece to adjust for the natural vision differential between eyes for viewing comfort. Camera-ready Trinocular models allow the user to view through both eyepieces while using the camera.
- Sturdy, 5-legged, easy-to-maneuver base slides in and under exam tables for optimal positioning
- Dual light intensity controls for both white and green LEDs
- Adjustable height center post with smooth suspension system
- S-Video and digital camera options available, trinocular version with built-in C-Mount Camera adapter

Complete with dust cover, power cord, manual, tool kit and 7 Year Limited Warranty, 1 year on electrical system (excludes bulbs and fuses). Optional accessories are available - please contact your Bovie Representative for additional information.

Item No.	Description
CS-205LED	Colpo-Master™ II Binocular Zoom Center Post Colposcope
CS-205LED-E	Colpo-Master™ II Binocular Zoom Center Post Colposcope, 220V
CS-205TLED	Colpo-Master™ II Camera-Ready Zoom Center Post Colposcope
CS-205TLED-E	Colpo-Master™ II Camera-Ready Zoom Center Post Colposcope, 220V

Colpo-Master™ is a registered trademark of United Products and Instruments, Inc; DBA UNICO



COLPO-MASTER™ III TRIPLE MAGNIFICATION COLPOSCOPES

Introducing the NEW Colpo-Master III Triple Magnification, LED Illuminated Center-Post Colposcopes. Select wheeled or pedestal base models.

Simple and easy to use, featuring brilliant, sharp optics and three magnification settings: 3.75x, 7.5x and 15x. Changing magnification is a simple rotation of the knob, while staying in focus!

The state of the art LED illumination provides bright, white light that will not burn out, or burn too hot!

Features:

- Three Magnification Settings: 3.75x, 7.5x and 15x
- Look straight through, no angled eyepieces
- Rheostat intensity controlled, bright LED illumination
- Paired 12.5x Widefield Eyepieces with rubber eye-guards
- 3-D Stereoscopic Viewing
- Dual Diopter, independent eye focusing controls
- Flip-down green filter, built into the head
- 300 mm Working Distance
- Extra-Wide 115 mm Interpupillary Distance Adjustment
- Center-Post stand with rolling, locking casters, or pedestal base
- Light weight, easy to maneuver
- Broad Range Universal Power Supply 90-240V
- Net Weight 23lbs, Shipping Weight 32lbs

Complete with dust cover, power cord, manual, tool kit and 7 Year Limited Warranty, 1 year on electrical system (excludes bulbs and fuses).

Item No.	Description
CS-305	Colpo-Master™ III Binocular Triple Magnification Colposcope with 5-Wheel Base
CS-305E	Colpo-Master™ III Binocular Triple Magnification Colposcope, 220V with 5-Wheel Base
CS-315	Colpo-Master™ III Binocular Triple Magnification Colposcope with Pedestal Base
CS-315E	Colpo-Master™ III Binocular Triple Magnification Colposcope, 220V with Pedestal Base
CS-8182	Digital Eyepiece Camera, 2.0 megapixel with Image Capture software for Windows XP/7

Colpo-Master™ is a registered trademark of United Products and Instruments, Inc; DBA UNICO