

Mobility and quality have now
become *one*. The new
HEINE DELTAone Dermatoscope.



30
YEARS
1st DERMATOSCOPE



Small on the outside. Big on the inside.

The new HEINE DELTAone Dermatoscope.



Extremely sharp image with the new achromatic HEINE optical system

Fast change from contact to non-contact examination

Fits into any pocket

Polarization eliminates reflections and glare

LEDHQ
LED NOW IN HEINE QUALITY.

Excellent color rendering due to LEDHQ

Digital documentation with the exclusive HEINE DERM app

- ∴ **Crystal clear image**, no matter if you prefer a closer or greater working distance
- ∴ **True field of view of 21 mm (lens 22mm) Ø**
- ∴ **Brightness control**: 3 levels of light intensity
- ∴ **Easily cleaned with spray disinfectant**
- ∴ **10x magnification**
- ∴ **No cross-contamination**: autoclavable contact plate with scale
- ∴ **Dioptré compensation**
- ∴ **Switch between polarized and non-polarized illumination mode at the touch of a button** = toggle function. For more information, please visit heine.com/en/Toggle

- ∴ **Rechargeable battery with Li-ion technology**: no memory effect; safely recharge the battery regardless of its state of charge
- ∴ **Small contact plate** for difficult-to-access pigmented lesions
- ∴ **2 possibilities of digital documentation**:
 - With the exclusive HEINE Derm App via iPhone* for model 6 or later and Android^{TM**}: Bodymap, patient folder, e-mail function and much more
 - Conventional digital images with many popular smartphones

DELTAone Standard Edition (black-silver) with contact plate with scale, USB cord with medical approved plug-in power supply, case	K-210.28.305
DELTAone White Edition (white-silver) with contact plate with scale, USB cord with medical approved plug-in power supply, case	K-211.28.305

∴ **Accessories and further articles under heine.com/en/DELTAone**

* The Apple iPhone is not included in the scope of delivery. Apple, the Apple logo and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc.

**Android is a trademark of Google LLC.

