

## ∴ HEINE® LoupeLight 2 with HR 2.5x Binocular Loupes

Perfect view for a precise diagnosis.

The combination of light and loupe delivers the ultimate in performance as high quality illumination maximises the image quality of the loupe. Featuring LED<sup>HQ</sup> illumination for homogeneous illumination of typ. 55,000 lux with accurate colour rendering that fills the entire field of view from edge-to-edge, the LoupeLight 2 is 25 % brighter than its predecessor. LED<sup>HQ</sup> – the new standard in LED illumination defining optimal light intensity, homogeneity and colour rendering for the most accurate diagnosis. Red is red, blue is blue.

**NEW!**

- ✓ Replaceable connecting cord
- ✓ Now 25 % brighter due to improved LED technology



**Replaceable connecting cord** with high-quality splashproof threaded connector for heavy duty use.

**LED<sup>HQ</sup>**

LED NOW IN HEINE QUALITY.

## ∴ HEINE® MicroLight 2

Ultralight, universal headlight.

For those who need a compact, lightweight illumination solution without magnification, we offer the MicroLight 2. Featuring LED<sup>HQ</sup> illumination for homogeneous illumination with true and natural colour rendering for every examination situation. Two mounting options offer individuality and easy handling. Illumination without compromises, in a lightweight design.



HEINE MicroLight 2  
on S-FRAME



HEINE MicroLight 2  
on Lightweight Headband

**NEW!**

- ✓ Replaceable connecting cord
- ✓ Now 25 % brighter due to improved LED technology



# HEINE® mPack mini

Accurate examinations made easy – all day long.

The new, high performance power source **HEINE mPack mini** for LoupeLight 2 and MicroLight 2 offers freedom from cables that tie the users to desk, wall or static power sources. With a weight of only 95g, this mobile power source is now approx. 70% smaller and lighter than its predecessor.

The pouch which is included in the delivery can be attached to the belt and the mPack mini can also be charged while clipped to the belt. Thanks to Li-ion battery technology the power source has an operating time of typ. 9 hours, in continuous operation at maximum power typ. 4 hours. The mPack mini can be charged flexibly via the USB plug, e.g. on a laptop/desktop or via USB charger.



## NEW!

- ✓ 70% smaller and lighter
- ✓ Operating time typ. 9 hours
- ✓ Charging time typ. 4 hours



Size comparison of HEINE mPack LL/mPack mini

## [ ORDERING INFORMATION ]

### LoupeLight 2 Sets with HR Binocular Loupes

HR 2.5x optics with i-View on S-FRAME, LoupeLight 2, detachable yellow filter, mPack mini, mPack mini Belt Clip, E4-USB plug-in transformer, protective lenses, 2 sterilisable levers, retaining cord, cleaning fluid, cleaning cloth, carrying case

340mm working distance	C-008.32.450
420mm working distance	C-008.32.451
520mm working distance	C-008.32.452

### MicroLight 2 Sets

MicroLight 2 on S-FRAME with mPack mini, mPack mini Belt Clip, E4-USB plug-in transformer, retaining cord and cleaning fluid in case

J-008.31.276

MicroLight 2 on Lightweight Headband with mPack mini, mPack mini Belt Clip, E4-USB plug-in transformer, in case

J-008.31.277

### mPack mini

mPack mini with Li-ion rechargeable battery and E4-USB plug-in transformer	X-007.99.650
mPack mini with Li-ion rechargeable battery, without E4-USB plug-in transformer	X-007.99.649
Connecting cord SC 1 (1.6 m) for LoupeLight 2 or MicroLight 2 / mPack mini	C-000.32.240

