

**Not All
Surgical Blades
Are Created Equal**



Trusted for Over 70 Years

featuring
Swann-Morton®


A Difference You Can Feel

The Cincinnati Surgical and Swann Morton Difference

A Different Feel


Our blades are manufactured with a unique cutting edge design. One side has a smooth razor edge that is supported by a micro serrated edge for a consistently sharp blade.

This grind delivers a tactile sensitivity that improves depth control while providing equal, smooth tissue margins. Many other blades have two razor edges. While initially very sharp, this edge can deteriorate faster.



Straight from the Manufacturer

- A commitment to quality
- Control over the entire manufacturing process
- Premium steel creates a uniform edge for a consistent cutting surface
- ISO certified for Quality, Energy and Environmental Management Systems
- 200% visual inspection resulting in zero defects
- In-house Sterilization
- Featured on the Discovery Channel's "How It's Made"



Carbon steel provides the best combination of initial sharpness & edge durability, but may stain with a solution wash.

Stainless steel contains a metal composition that keeps it from staining when used with a solution. Known for its lateral pressure durability.

Dental Blades

Standard and Specialty - Sterile, packaged 100 per box



Swann Morton
00SM11: carbon
01SM11: stainless
Cincinnati Surgical
0011: carbon
0111: stainless



Swann Morton
00SM11P: carbon
01SM11P: stainless



Swann Morton
00SM12: carbon
01SM12: stainless
Cincinnati Surgical
0012: carbon
0112: stainless



Swann Morton
00SM12D*: carbon
01SM12D*: stainless
Cincinnati Surgical
0112D*: stainless
*same as 12B



Swann Morton
00SM15: carbon
01SM15: stainless



Swann Morton
00SM15C: carbon
01SM15C: stainless
Cincinnati Surgical
0115C: stainless



Swann Morton
00SM15T: carbon
01SM15T: stainless
Cincinnati Surgical
0015: carbon
0115: stainless

Miniature- Sterile, packaged 25 per box



Swann Morton
01SM64: stainless



Swann Morton
01SM65: stainless



Swann Morton
01SM67: stainless



Swann Morton
01SM69: stainless

Disposable Scalpels- Sterile, stainless, packaged 10 per box



Size 11: SM0503
Size 12: SM0504
Size 12D: SM0518



Size 15: SM0505
Size 15T: SM0592

Size 11: 0511
Size 12: 0512
Size 15: 0515

About the Manufacturer

Since 1938, Cincinnati Surgical has been recognized as an industry leader in manufacturing high-quality surgical blades, scalpels, and safety products. Our manufacturing heritage in the United States (1938-1961) and Sheffield, England (1962-Present) has now come under one roof in Sheffield through a historic joint manufacturing agreement with Swann Morton, signed in 2015. By combining our manufacturing capabilities, we've expanded our family of brands and brought the global leader in surgical blade manufacturing to the US market.

From the quality of the steel to our ISO manufacturing standards and visual inspection of each blade, we are able to provide a consistent, zero defect product line that will ensure durability and consistency each time it is used. With our extensive range of products and ability to ship same day, we offer our customers a product they know they can trust when quality counts.

Our number one focus is to produce the highest quality blade available. Because surgical blades are a Class 1 Medical Device they do not require FDA inspections, whether imported or repackaged. This can lead to diminished quality control by some in the industry. But don't worry because "We've Got You Covered" through every step of production. From every batch of steel meeting our strict specifications all the way through to our in-house sterilization process.

To **feel the difference** for yourself purchase or sample these blades by contacting:



Trusted for Over 70 Years

featuring
Swann-Morton®

www.cincinnati-surgical.com

csblades@cincinnati-surgical.com

800-544-3100