

# Not All Blades Are Created Equal



featuring:  
**Swann-Morton®**  
CincinnatiSurgical.com

## A Different Feel

Our blades are manufactured with a unique cutting edge design. One side has a smooth razor edge that is supported by a micro serrated edge for a consistently sharp blade.

This grind delivers a tactile sensitivity that improves depth control while providing equal, smooth tissue margins. Many other blades have two razor edges. While initially very sharp, this edge can deteriorate faster.



Since 1938, Cincinnati Surgical has been recognized as an industry leader in manufacturing high-quality surgical blades, scalpels, and safety products. Our manufacturing heritage in the United States (1938-1961) and Sheffield, England (1962-Present) has now come under one roof in Sheffield through a historic joint manufacturing agreement with Swann Morton, signed in 2015. By combining our manufacturing capabilities, we've expanded our family of brands and brought the global leader in surgical blade manufacturing to the US market.

From the quality of the steel to our ISO manufacturing standards and visual inspection of each blade, we are able to provide a consistent, zero defect product line that will ensure durability and consistency each time it is used. With our extensive range of products and ability to ship same day, we offer our customers a product they know they can trust when quality counts.


Our number one focus is to produce the highest quality blade available. Because surgical blades are a Class 1 Medical Device they do not require FDA inspections, whether imported or repackaged. This can lead to diminished quality control by some in the industry. But don't worry because "We've Got You Covered" through every step of blade production. From every batch of steel meeting our strict specifications all the way through to our in-house sterilization process.

To **feel the difference** for yourself purchase or sample these blades by contacting:

Tracy Sterling  
tracy@cincinnatiSurgical.com  
513-520-0542

## Straight from the Manufacturer

- A commitment to quality
- Control over the entire manufacturing process
- Premium steel creates a uniform edge for a consistent cutting surface
- ISO certified for Quality, Energy and Environmental Management Systems
- 200% visual inspection resulting in zero defects
- In-house Sterilization
- Featured on the Discovery Channel's "How It's Made"



Carbon steel provides the best combination of initial sharpness & edge durability, but may stain with a solution wash.

Stainless steel contains a metal composition that keeps it from staining when used with a solution. Known for its lateral pressure durability.

# Veterinary Blades

**Standard and Specialty** - Sterile, packaged 100 per box



**Swann Morton**  
 00SM11: carbon  
 01SM11: stainless  
**Cincinnati Surgical**  
 0011: carbon  
 0111: stainless



**Swann Morton**  
 00SM11P: carbon  
 01SM11P: stainless



**Swann Morton**  
 00SM12: carbon  
 01SM12: stainless  
**Cincinnati Surgical**  
 0012: carbon  
 0112: stainless



**Swann Morton**  
 00SM12D\*: carbon  
 01SM12D\*: stainless  
**Cincinnati Surgical**  
 0112D\*: stainless  
 \*same as 12B



**Swann Morton**  
 00SM15C: carbon  
 01SM15C: stainless  
**Cincinnati Surgical**  
 0115C: stainless



**Swann Morton**  
 00SM15T: carbon  
 01SM15T: stainless  
**Cincinnati Surgical**  
 0015: carbon  
 0115: stainless



**Swann Morton**  
 00SM20: carbon  
 01SM20: stainless  
**Cincinnati Surgical**  
 0120 stainless



**Swann Morton**  
 00SM21: carbon  
 01SM21: stainless  
**Cincinnati Surgical**  
 0121: stainless



**Swann Morton**  
 00SM22: carbon  
 01SM22: stainless  
**Cincinnati Surgical**  
 0022 carbon  
 0122 stainless

**Miniature** - Sterile, packaged 25 per box



**Swann Morton**  
 01SM64: stainless



**Swann Morton**  
 01SM65: stainless



**Swann Morton**  
 01SM67: stainless



**Swann Morton**  
 01SM69: stainless

**Disposable Scalpels**- Sterile, stainless, packaged 10 per box



Size 11: SM0503  
 Size 12: SM0504  
 Size 12D: SM0518  
 Size 15: SM0505  
 Size 20: SM0506  
 Size 21: SM0507  
 Size 22: SM0508



Size 11: 0511  
 Size 12: 0512  
 Size 15: 0515  
 Size 20: 0520  
 Size 21: 0521  
 Size 22: 0522