



COVIDIEN

positive results for life™


Maxon™

Monofilament Synthetic Absorbable Suture

**THE STRENGTH AND FLEXIBILITY YOU REQUIRE
YOUR CHOICE FOR LONG-LASTING
WOUND SUPPORT**

- Long-term monofilament absorbable suture
- Unequalled strength during the critical stages of healing*
- Six weeks of wound support
- Monofilament construction provides smooth passage through tissue, minimizing tissue drag and trauma
- Maxon™ looped sutures provide double the strength of a single strand and can reduce closure time by 50%
- Now at least 36% more pliable than “old” Maxon™ on sizes 2-0, 0 and 1

**Maxon™ Suture
36% Pliability Improvement (Size 0)**



Frequent Uses

General soft tissue approximation including fascia closure, Orthopedic surgery, Pediatric Cardiovascular surgery and patients with compromised wound healing.

Indications


Maxon™ and Maxon™ CV synthetic absorbable sutures are indicated for use in general soft tissue approximation and/or ligation, including use in pediatric cardiovascular tissue, where growth is expected to occur, and in peripheral vascular surgery. Maxon™ and Maxon™ CV synthetic absorbable sutures are not indicated for use in adult cardiovascular tissue, ophthalmic surgery, microsurgery and neural tissues.

Specifications†


Composition	Glycolide and Trimethylene Carbonate
Sizes	Available in USP sizes to 7-0 to 1
Color	Green and Clear
Absorption Profile	180 Days
Tensile Strength	80% of initial tensile strength retained at 1 week, 75% at 2 weeks, 65% at 3 weeks, 50% at 4 weeks and 25% at 6 weeks post implant

Maxon™ Is Stronger Than The Competition During The Critical Wound Healing Period

Maxon™ Suture – Initial Strength (Size 0)



*Strength During Critical Wound Healing (Size 0)



Learn more: www.covidien.com/syneture
or place your order at: 1 800 722 8772.



COVIDIEN, COVIDIEN with logo and "positive results for life" are trademarks of Covidien AG.
© 2008 Covidien AG or its affiliates. All rights reserved.
378969 SS 2.5M 08.08 CSS000170i

†Refer to the package insert for complete product details