



PENCILS


The standard Bovie pencils feature button, rocker, or foot control, available with and without holster and scratch pad. All disposable pencils are packaged individually sterile and are offered with stainless steel or Resistick II coated electrodes. Additionally, Bovie offers a reusable pencil (ESPR). This pencil is validated to forty (40) autoclave cycles and is available sterile or non-sterile.

	ESPR (Top 40)* autoclavable reusable pencil, non-sterile - 1/box (*electrode not included) <i>Pencils use a standard 3-prong connector</i> ESPRS (Top 40) autoclavable reusable pencil, sterile - 5/box <i>Pencils use a standard 3-prong connector</i>		 3-prong connector
	ESP1 disposable pencil, sterile	CONFIGURATIONS ESP1: pencil, blade - 50/box ESP1H: pencil, blade, holster - 40/box ESP1HS: pencil, blade, holster, scratch pad - 40/box ESP1N: pencil, needle - 50/box ESP1HN: pencil, needle, holster - 40/box <i>Pencils use a standard 3-prong connector</i>	
	ESP6 disposable rocker switch pencil, sterile	CONFIGURATIONS ESP6: pencil, blade - 50/box ESP6H: pencil, blade, holster - 40/box ESP6HS: pencil, blade, holster, scratch pad - 40/box ESP6N: pencil, needle - 50/box ESP6HN: pencil, needle, holster - 40/box <i>Pencils use a standard 3-prong connector</i>	
	ESP7 disposable foot-control pencil, sterile	CONFIGURATIONS ESP7: pencil, blade - 50/box ESP7H: pencil, blade, holster - 40/box ESP7HS: pencil, blade, holster, scratch pad - 40/box ESP7N: pencil, needle - 50/box ESP7HN: pencil, needle, holster - 40/box <i>Pencils use a single pin connector</i>	


RESISTICK II™ DISPOSABLE PENCILS			
ESP1T	ESP1TN	ESP6T	ESP6TN
ESP1T: disposable pencil, ES01T blade - 50/box		ESP6T: disposable pencil, ES01T blade - 50/box	
ESP1TN: disposable pencil, ES02T needle - 50/box		ESP6TN: disposable pencil, ES02T needle - 50/box	

PENCILS & ACCESSORIES


Continued...




A901
autoclavable reusable 3-button pencil for
Aaron 900™, A940™ & 950™ - comes with
generator or 1/box



A910
disposable handpiece sheath, non-sterile
- 100/box



ESPH
disposable pencil holster, sterile - 40/box



ESSP
disposable scratch pad - 40/box


SUCTION COAGULATORS

Continuing to expand its electrosurgical offering, Bovie Medical Corporation is pleased to introduce six Suction Coagulators. They are compatible with Bovie® and Aaron® generators, as well as, those of most other companies. Featuring three different diameters, 6" flexible shaft and hand or foot activation, these suction coagulators allow fluid evacuation with or without simultaneous coagulation. These single use devices are individually packaged sterile, ten per shelf pack.



HANDSWITCH SUCTION COAGULATORS


Catalog Number	Description	Tube Diameter (French Size)	Cable Length
SCH08	Disposable	08 Fr	3m
SCH10	Disposable	10 Fr	3m
SCH12	Disposable	12 Fr	3m



	12	10	8
French Gauge	12	10	8
Diameter (mm)	4	3.3	2.7
Diameter (inches)	0.158	0.131	0.105

FOOTSWITCH SUCTION COAGULATORS

Catalog Number	Description	Tube Diameter (French Size)	Cable Length
SCF08	Disposable	08 Fr	3m
SCF10	Disposable	10 Fr	3m
SCF12	Disposable	12 Fr	3m



STANDS AND WALL MOUNT			
	A812-C mobile stand for Aaron® 800, 800-EU, 900, 940 & 950 (includes bottom tray, top tray, clamp, and 12" height extension) 47" high		ESMS-C mobile stand for Aaron® 1250, 2250 & 3250 (includes stand & bottom tray) 24" wide x 24" deep x 31" high
	A812 mobile stand for Aaron® 800, 800-EU, 900, 940 & 950, 35" high		ESMS mobile stand for Aaron® 1250, 2250 & 3250 24" wide x 24" deep x 31" high
	A812-BT bottom tray for A812 20" wide x 11" deep		ESMS-B bottom tray for ESMS 12" wide x 7.75" deep
	A808-T top tray and clamp for A812 13.5" wide x 9.5" deep		BV-IDS-CS mobile stand for Aaron® 1250, 2250 & 3250 16" wide x 20" deep x 35" high
	A812-E 12" height extension for A812		BV-ICON-CS* mobile stand for Icon Gi generator 18" wide x 18" deep x 34" high (without generator) *Icon Gi generator sold separately—see page 11
	A837 wall mount kit for Aaron® 800-EU, 800, 900, 940 & 950		A813 table-top stainless steel stand for Aaron® 800, 800-EU, 900, 940, 950 & 1200

STANDS, FOOTSWITCHES & WALL MOUNT KIT




FOOTSWITCHES	
	A803 Aaron® 800, 800-EU & 900 and 940 footswitch
	A1203 and A1203W (Waterproof) Aaron® 950 & 1200 footswitch
	BV-1253BGi Icon Gi footswitch
	BV-1253B Aaron® 1250, 2250 & 3250 footswitch
	BV-1254B Aaron® 2250 & 3250 bipolar footswitch



DISPERSIVE ELECTRODES, CORDS, ADAPTERS

Aaron disposable dispersive electrodes with super adhesive gel are designed for single use with safety and quality built in. Where reusability is requested, a reusable electrode and cable are available.

DISPOSABLE GROUNDING PADS AND CORDS	
	A1202 or ESRS solid dispersive electrode A1202 - 5/pack ESRS - 50/box
	ESRSC solid adult return electrode w/2.8 m cable 50/box
	ESRE split adult return electrode 50/box
	ESREC split adult return electrode w/2.8 m cable - 50/box

ADAPTERS	
	A905A (shown) and A1205A adapter for connecting foot- switch pencil to the Aaron® 900, 940 and 950
	A905B adapter for 1/32" (0.8mm) electrode -fits most ESU pencils
	A905B-1 adapter for 1/20" (1.25mm) electrode -fits most ESU pencils
	A905EL adapter for 1/16" electrode into 3/32" collet
	A1255A adapter plug for connecting footswitching pencil to the Aaron® 1250, 2250, 3250 and most ESUs

REUSABLE GROUNDING PADS AND CORDS	
	A802EU metal dispersive electrode for Aaron® 800-EU & 900 and 940 - 1/pouch
	A1204 metal dispersive electrode and cord for Aaron® 950 - 1/pouch
	A1204P replacement metal dispersive electrode for use with A1204C and A1254C - 1/pouch
	A1202C connecting cord for A1202 or ESRS dispersive electrode for Aaron® 950 - 1/pouch
	A1252C connecting cord for ESRE and ESRS to Aaron® 1250, 2250, 3250 and most ESUs - 1/pouch
	A1204C replacement cord for Aaron® 950 - 1/pouch
	A1254C same as A1204C but for Aaron® 1250, 2250 and 3250 - 1/pouch