




DISPOSABLE ELECTRODES

All disposable electrodes utilize the standard 3/32" (2.3 mm) stainless steel shaft. All molded electrodes feature safety grip insulators combining patient and user safety with easy insertion into and removal from the surgical pencil. disposable electrodes are manufactured to the highest standards. They come individually packaged sterile*.



BLADES		BALLS	NEEDLES
ES01 standard blade 2.75" (6.985cm) 0.75" (1.90cm) - 50/box ES55 extended blade 4" (10.16 cm) - 25/box ES04 extended blade 6" (15.24cm) - 25/box ES37 modified blade 2.75" (6.985cm) - 25/box ES54 extended modified blade 4" (10.16cm) - 25/box ES39 extended modified blade 6" (15.24cm) - 25/box ES18 angled blade 2.5" (6.35cm) - 5/box	A804* sharp dermal tip for use with Aaron® 800, 900, 940, & 950 - 100/box A805 sharp dermal tip for use with Aaron® 800, 900, 940, & 950, sterile - 50/box A806* blunt dermal tip for use with Aaron® 800, 900, 940, & 950 - 100/box A807 blunt dermal tip for use with Aaron® 800, 900, 940, & 950, sterile - 50/box	ES06 extended ball 0.118" (3mm) - 5/box ES07 extended ball 0.197" (5mm) - 5/box ES50 1/8" ball 0.125" (3.175mm) - 25/box ES20 3/16" ball 0.188" (4.775mm) - 25/box ES21 5/32" ball 0.15625 (3.962mm) - 25/box	ES02 standard needle 2.75" (6.985cm) - 25/box ES03 extended needle 6" (15.24cm) - 25/box ES38 modified needle 2.75" (6.985cm) - 25/box ES40 extended modified needle 6" (15.24cm) - 25/box

* A804 & A806 are non-sterile and not packaged individually.

REUSABLE ELECTRODES



All disposable electrodes utilize the standard 3/32" (2.3 mm) stainless steel shaft. All molded electrodes feature safety grip insulators combining patient and user safety with easy insertion into and removal from the surgical pencil. Reusable electrodes can be steam autoclaved up to 25 times. Reusable loop electrodes can be steam autoclaved up to five times (see page 20).





ELECTRODES	BLADES	BALLS	NEEDLES
A811 reusable blunt electrode - 1/box A830 angled sharp electrode - 1/box A836 extended angled sharp electrode - 1/box	ES01R standard blade - 1/box ES04R extended blade - 1/box ES18R angled blade - 1/box	A831 angled ball - 1/box A832M 3/32" ball - 1/box ES06R extended 3 mm ball - 1/box ES07R extended 5 mm ball - 1/box ES20R 3/16" ball - 1/box ES21R 5/32" ball - 1/box	A833 straight needle - 1/box A834 angled fine (epilation) needle - 1/box A834T angled Tungsten fine needle - 1/box ES02R standard needle - 1/box ES03R extended needle - 1/box






RESISTICK II™ COATED ELECTRODES


Resistick II™ coated blade, ball and needle electrodes are expertly coated with PTFE (polytetrafluoroethylene). They move more effectively through tissue because they resist eschar build-up making them much more efficient during surgical procedures. Just use a damp sponge or cloth to almost effortlessly clean the debris off the Resistick II™ electrodes.



	ES50T
Standard 3mm Ball Electrode 2" (5.08cm) - 12/box	
	ES21T
Standard 4mm Ball Electrode 2" (5.08cm) - 12/box	
	ES20T
Standard 5mm Ball Electrode 2" (5.08cm) - 12/box	
	ES06T
Extended 3mm Ball Electrode 5" (12.70) - 12/box	
	ES07T
Extended 5mm Ball Electrode 5" (12.70cm) - 12/box	

	SS01T
Super Sharp Modified Coated 2" (5.08cm) Needle Electrode 12/bx	
	SS02T
Super Sharp Modified Coated 2.5" (6.35cm) Needle Electrode 12/bx	
	ES02T
Standard Needle Electrode 2.75" (6.98cm) - 12/box	
	ES38T
Modified Needle Electrode 2.75" (6.98cm) - 12/box	
	ES56T
Extended Needle Electrode 4" (10.16cm) - 12/box	
	ES57T
Extended Modified Needle Electrode 4" (10.16cm) - 12/box	
	ES03T
Extended Needle Electrode 6" (15.24cm) - 12/box	
	ES40T
Extended Modified Needle Electrode 6" (15.24cm) - 12/box	

	ES58T
Blade Electrode 2.5" (6.35cm) - 12/box	
	ES59T
Modified Blade Electrode 2.5" (6.35cm) - 12/box	
	ES01T
Standard Blade Electrode 2.75" (6.98cm) - 12/box	
	ES37T
Modified Blade Electrode 2.75" (6.98cm) - 12/box	
	ES18T
Angled Blade 45° Electrode 2.75" (6.98cm) - 12/box	
	ES55T
Extended Blade Electrode 4" (10.16cm) - 12/box	
	ES54T
Extended Modified Blade Electrode 4" (10.16cm) - 12/box	
	ES04T
Extended Blade Electrode 6" (15.24cm) - 12/box	
	ES39T
Extended Modified Blade Electrode 6" (15.24cm) - 12/box	



DISPOSABLE LOOP ELECTRODES

All Aaron® LLETZ Loop* and Square electrodes* feature tungsten wire for superior shape and integrity throughout the excision procedure. They use a standard $\frac{3}{32}$ " (2.3 mm) shaft. All molded loops are packaged in a patented protective tray to prevent damage during shipping.





LOOPS

ES08 5 mm X 5 mm loop - 5/box	ES31 20 mm X 20 mm loop - 5/box	ES47 25 mm x 10 mm loop - 5/box
ES09 10 mm X 10 mm loop - 5/box	ES42 20 mm X 12 mm loop - 5/box	ES49 20 mm X 10 mm loop - 5/box
ES10 15 mm X 8 mm loop - 5/box	ES43 15 mm x 15 mm loop - 5/box	ES51 13 mm X 8 mm loop - 5/box
ES11 15 mm X 10 mm loop - 5/box	ES44 15 mm x 5 mm loop - 5/box	ES52 10 mm X 8 mm loop - 5/box
ES12 20 mm X 8 mm loop - 5/box	ES45 13 mm x 13 mm loop - 5/box	ES53 15 mm X 12 mm loop - 5/box
ES13 20 mm X 15 mm loop - 5/box	ES46 10 mm x 5 mm loop - 5/box	

SQUARES



ES41 5 mm X 10 mm square - 5/box	ES14 5 mm X 5 mm square - 5/box	ES15 10 mm X 4 mm square - 5/box	ES16 10 mm X 8 mm square - 5/box
--	---	--	--






REUSABLE LOOPS/SQUARES, ELECTRODES

All Aaron® LLETZ Loop* and Square* electrodes feature tungsten wire for superior shape and integrity throughout the excision procedure. They use a standard 3/32" (2.3 mm) shaft. All molded loops are packaged in a patented protective tray to prevent damage during shipping.



LOOPS		
ES08R 5 mm X 5 mm loop - 1/box	ES31R 20 mm X 20 mm loop - 1/box	ES47R 25 mm x 10 mm loop - 1/box
ES09R 10 mm X 10 mm loop - 1/box	ES42R 20 mm X 12 mm loop - 1/box	ES49R 20 mm X 10 mm loop - 1/box
ES10R 15 mm X 8 mm loop - 1/box	ES43R 15 mm x 15 mm loop - 1/box	ES51R 13 mm X 8 mm loop - 1/box
ES11R 15 mm X 10 mm loop - 1/box	ES44R 15 mm x 5 mm loop - 1/box	ES52R 10 mm X 8 mm loop - 1/box
ES12R 20 mm X 8 mm loop - 1/box	ES45R 13 mm x 13 mm loop - 1/box	ES53R 15 mm X 12 mm loop - 1/box
ES13R 20 mm X 15 mm loop - 1/box	ES46R 10 mm x 5 mm loop - 1/box	

SQUARES			
ES41R 5 mm X 10 mm square - 1/box	ES14R 5 mm X 5 mm square - 1/box	ES15R 10 mm X 4 mm square - 1/box	ES16R 10 mm X 8 mm square - 1/box



DISPOSABLE ARTHROSCOPIC ELECTRODES

AR03 – Menisectomy, sterile - 5/box
AR02 – Acromioplasty, sterile - 5/box

AR00 – hook 90, sterile - 5/box
AR01 – hook 45, sterile - 5/box




REUSABLE SUPERCUT™ TUNGSTEN NEEDLES

ES60R – needle super fine 3 cm, non-sterile - 1/box
ES62R – needle super fine 4.5 cm, non-sterile - 1/box

DISPOSABLE SUPERCUT™ TUNGSTEN NEEDLES

ES60 – needle super fine 3 cm, sterile - 5/box
ES61 – modified needle super fine 3 cm, sterile - 5/box
ES62 – needle super fine 4.5 cm, sterile - 5/box
ES63 – modified needle super fine 4.5 cm, sterile - 5/box

TUNGSTEN LOOPS, NEEDLES & ARTHROSCOPIC ELECTRODES



DISPOSABLE SHORT TUNGSTEN .008" WIRE LOOPS

ES22-8 – 1/4" (6.35mm) short shaft loop, sterile - 5/box
ES23-8 – 5/16" (7.94mm) short shaft loop, sterile - 5/box
ES24-8 – 3/8" (9.53mm) short shaft loop, sterile - 5/box
ES25-8 – 1/2" (12.7mm) short shaft loop, sterile - 5/box
ES26-8 – 3/4" (19.05mm) short shaft loop, sterile - 5/box

DISPOSABLE SHORT TUNGSTEN .015" WIRE LOOPS

ES22 – 1/4" (6.35mm) short shaft loop, sterile - 5/box
ES23 – 5/16" (7.94mm) short shaft loop, sterile - 5/box
ES24 – 3/8" (9.53mm) short shaft loop, sterile - 5/box
ES25 – 1/2" (12.7mm) short shaft loop, sterile - 5/box
ES26 – 3/4" (19.05mm) short shaft loop, sterile - 5/box

DISPOSABLE SUPERCUT™ TUNGSTEN NEEDLES

ES61HS – modified fine needle with Heatshrink, sterile - 5/box



REUSABLE SHORT TUNGSTEN .015" WIRE LOOPS

ES22R – 1/4" (6.35mm) short shaft loop, non-sterile - 1/box
ES23R – 5/16" (7.94mm) short shaft loop, non-sterile - 1/box
ES24R – 3/8" (9.53mm) short shaft loop, non-sterile - 1/box
ES25R – 1/2" (12.7mm) short shaft loop, non-sterile - 1/box
ES26R – 3/4" (19.05mm) short shaft loop, non-sterile - 1/box



REUSABLE BIPOLAR FORCEPS

Reusable Bipolar Forceps are available in micro (.5 mm), 1.0 mm, and 2.0 mm tip size and in lengths ranging from 3.5" to 8". These forceps are guaranteed for twenty autoclave cycles. A reusable bipolar cord with quality silicone construction is also available.

		A820 7.25" bayonet, fine (1.0 mm), smooth - 1/box	A821 8" bayonet, medium (2.0 mm), smooth - 1/box	A822 7" straight, fine (1.0 mm), smooth - 1/box
A823 7" curved, fine (1.0 mm), smooth - 1/box	A824 4.5" curved, micro (.5 mm) - 1/box	A825 4.75" straight, medium (2.0 mm), smooth - 1/box	A826 4.5" straight, micro (.5 mm), smooth - 1/box	A840 4" Tenzel, micro (.5 mm), coated - 1/box
A841 3.5" McPherson, curved, micro (.5 mm), coated - 1/box	A842 3.5" McPherson, straight, micro (.5 mm), coated - 1/box	A843 4" Tenzel, micro (.5 mm), uncoated - 1/box	A844 3.5" McPherson, curved, micro (.5 mm), uncoated - 1/box	A845 3.5" McPherson, straight, micro (.5 mm), uncoated - 1/box
A827 forceps cord - 1/pouch fits Aaron® 800, 1250, 2250, 3250 & most ESUs	A827BP forceps cord - 1/pouch fits Aaron® 900, 940 & 950	A827EU forceps cord - 1/pouch fits Aaron® 800-EU, 900, 940 & 950	A827F forceps cord - 1/pouch fits Aaron® 1250(only)	A827V forceps cord - 1/pouch fits Aaron® 2250 & 3250