


FLEXIBLE LIGHTS AND STYLETS

Bovie manufactures two groups of flexible lighting products. The first is available both sterile and non-sterile in three different sizes 5", 10", and 15" (refers to the flexible shaft portion of the light). These lights are for use in procedures where additional light is needed. The second group of flexible lights is our lighted stylets, for difficult intubations.



FLEXIBLE LIGHTS		LIGHTED STYLETS	
NS05 * Surch-Lite™ 5" non-sterile - 1 ea	ST05 Surch-Lite™ 5" sterile - 10/box	SLOT oro-tracheal stylet, sterile - 1/box	PDOT pediatric oro-tracheal stylet, sterile - 1/box
NS10 * Surch-Lite™ 10" non-sterile - 1 ea	ST10 Surch-Lite™ 10" sterile - 3/box	NLOT * oro-tracheal stylet, non-sterile - 1/box	PDON * pediatric oro-tracheal stylet, non-sterile - 1/box
NS15 Surch-Lite™ 15" non-sterile - 1 ea	ST15* Surch-Lite™ 15" sterile - 3/box		

UV WOODS LIGHT, SPECIALTY LIGHTS, PENLIGHTS



WOODS LIGHT

UV59

hand-held woods light, excellent for detecting certain skin disorders and rape detection procedures - 1/box



SPECIALTY LIGHTS



6800

tech light - 1/box




6801

white slit - 1/box



6802

Fluoro-Dot™ - 1/box




6803

Fluoro-Slit™ - 1/box



6804

1mm light - 1/box



INDIVIDUAL PENLIGHTS

1135

doctor's replaceable
penlight - 1/box



1125

doctor's disposable
penlight - 1/box



PENLIGHT PACKS

6303

cobalt penlight - 3/pack



666R

doctor's replaceable battery penlight
- 6/pack



6666

doctor's disposable penlight - 6/pack



66RN

penlight w/pupil gauge - 6/pack



7156

cobalt filter only - 1/pack

