



RELIABLE
Physician's
Chair Scales
for
Patients
with Special
Health Care
Needs



Footrest lifts out of the way for safe seating and convenient transport of scale while the lift-away arms allow patients side or front access. Oversized wheels facilitate movement of scale. Complete with heavy-duty understructure, these scales come fully assembled and ready to use. Black vinyl upholstery is complemented by the metal support structure.

MECHANICAL CHAIR SCALE

Mechanical chair scales are available in pounds or pounds/kilo models. Die-cast, easy-to-read beams are readable from front and back.

Capacity - 475: 400 lb x 4 oz

4751: 400 lb x 4 oz/180 kg x 100 g

DIGITAL CHAIR SCALE

Digital chair scales are available in pound or kilo models. Easy-to-read digital display rapidly shows weight to an accuracy of .2 lb (.1 kg). Battery power enhances portability.

Capacity - 6475: 400 lb x .2 lb

6475KG: 180 kg x .1 kg

Display: .75" high-contrast LCD

Power: 6 "AA" size Alkaline batteries

(not included)

Option: 115 VAC power adapter



Cardinal Scale reserves the right to improve, modify, or enhance features and specifications without prior notice.

PH: 417-673-4631 • FAX: 417-673-5001 P.O. BOX 151, WEBB CITY, MO 64870 USA

WEB SITE: http://www.detectoscale.com