MetriCide™ OPA Plus

MetriCide OPA Plus is the ideal choice for healthcare professionals who are concerned about the cost and effectiveness of high-level disinfectants used in reprocessing endoscopes and other heat sensitive semi-critical medical devices.



Features & Benefits:

- Manual reprocessing: can be reused for up to 30 days in 12 minutes at 20°C No other OPA lasts longer manually
- Automated reprocessing (AER): can be reused for up to 14 days in 5 minutes at 25°C
- Surfactant free formula won't cloud scope lenses

Efficacy

The following table indicates the spectrum of activity as demonstrated by testing of MetriCide OPA Plus using prescribed test methods. Bactericidal, Fungicidal, Tuberculocidal Tests all passed in 5 minutes at 25°C.

Microorganism **Viruses VEGETATIVE ORGANISMS NON-ENVELOPED ENVELOPED** Staphylococcus aureus Poliovirus Type 1 Influenza Virus Type A [Hong Kong] Salmonella choleraesuis Rhinovirus Type 42 HIV-1 Pseudomonas aeruginosa Hepatitis A (AER) Herpes simplex Type 1 Mycobacterium tuberculosis var: bovis Duck Hepatitis B Avian Influenza A

FUNGI

Trichophyton mentagrophytes

Ordering Information

Part No.	Description	Packaging
10-6000	MetriCide OPA Plus Solution 1 gallon	4 bottles/case
10-602	MetriCide OPA Plus Solution Test Strips - 100 Strips per bottle	2 bottles/case

For more information on MetriCide™ OPA Plus, visit Metrex.com or call 800-841-1428.

