



The Schiller 3 Liter Calibration Syringe is designed to simplify the process of spirometry calibration and ensure that your spirometer devices are functioning correctly (ATS2005 Guidelines).



GENERAL SPECIFICATIONS

2.100037L

Calibration Syringe Universal Adapter

3 Liters \pm 0,5%

Weight: 3,8 lbs.

Dimensions: 23" x 4,1"

Mouthpiece dimension: 1,25"



Options

4.430181

Calibration Syringe Universal Adapter



OPERATION

Schiller recommends yearly syringes calibration in order to ensure accuracy and conform to the ATS2005 Standardization of Spirometry quality control guidelines

SCHILLER CALIBRATION SYRINGE