III Hillrom...





Welch Allyn® Connex® Integrated Wall System

The information and tools you need, all in one place.

The Welch Allyn[®] Connex[®] Integrated Wall System is an advanced vitals management and physical assessment solution that combines the exam tools you need in one convenient system. Take and electronically capture accurate patient vital signs and perform basic diagnostic exams, all while saving the valuable time and space.

Studies show automated blood pressure averaging can provide additional decision support to improve the accuracy of hypertension diagnoses.*

- Capture numerous readings in succession and average them automatically for a more complete and accurate picture of your patient's heart health
- Automated office blood pressure devices can be programmed to record and average multiple blood pressure readings — with or without a clinician present — and can provide a more accurate blood pressure measurement than manual auscultation**

*Robert Smith, MD. Blood pressure averaging methodology: Decreasing the risk of misdiagnosing hypertension

**5 Muntner, Paul et al. "Measurement of Blood Pressure in Humans: A Scientific Statement From the American Heart Association." Hypertension 73, no. 5. pg.e35. (May 1, 2019). Doi:10.1161/ hyp.000000000000087



Configure your wall system based on your workflow:



Integrated otoscope and ophthalmoscope (choose from a variety of 3.5 V options)



Adult, pediatric and neonatal blood pressure modes featuring SureBP[®] technology (15-second NIBP and pulse rate)



Choice of Masimo® or Nellcor® pulse oximetry



Immediately send accurate patient vitals directly to your EMR.

- Reduce manual transcription errors in connected facilitiested
- Help improve patient safety and reduce the time spent documenting vitals
- View patient data anytime and anywhere on your network

Expand and upgrade your system with a flexible design.

- Add new parameters and devices as your needs change and new technology becomes available
- Protect your investment and scale your solution to your care environment's needs
- Perform software-driven repairs and upgrades

Enhance clinical decision-making, improve patient safety and save valuable time with integrated vital signs management and physical assessment tools.

- Configure your wall system based on your workflow
- Eliminate frustrating and time-consuming searches for individual pieces of equipment
- Prevent tools from being lost, stolen or misplaced
- Improve information visibility for your patients to enable easier communication and promote compliance
- Save valuable floor space and minimize device wear and tear
- Get the accurate data you need to enable data-driven decision-making for improved patient car



Integrated thermometry and probe cover holder for the SureTemp[®] Plus oral thermometer and/or Braun ThermoScan[®] PRO 6000 ear thermometry (choose one or both)

Large, easy-to-read color touchscreen

- Show patients their results as they happen for improved communication and patient compliance
- NIBP interval monitoring is available with customizable interval programs
- Manually enterl measurements like weight, height, respiration rate, pain — and document modifiers like body position and O₂ therapy

 Functions are grouped to easily comply with Meaningful Use requirements


Wall-mounted cord/cuff management system reduces tangled, hanging cords and provides organized storage for extra blood pressure cuffs

Hillrom.

Let Hillrom help minimize interruptions to patient care and unexpected costs caused by equipment that is out for repair or maintenance. Our SmartCare[™] Protection and Protection+ programs deliver the right mix of services to help you reduce equipment downtime and keep devices in the hands of your caregivers.

- Experience fewer interruptions in care
- Maximize the life of your devices and systems
- Reduce repair costs
- Stay compliant with local and national accreditations

Contact your Hillrom representative for more details on how our proactive services can provide you with peace of mind.



For more information, please contact your local Hillrom representative or visit us at hillrom.com.

hillrom.com

4341 State Street Rd., Skaneateles Falls, NY 13153

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

Braun and ThermoScan are registered trademarks of Braun, GmbH. Masimo is a registered trademark of Masimo Corporation. Nellcor and OxiMax are registered trademarks of Nellcor Puritan Bennett, Inc.

© 2021 Welch Allyn, Inc. ALL RIGHTS RESERVED. APR124001 Rev 1 16-AUG-2021 ENG-US