



## Welch Allyn KleenSpec® Single Use LED Vaginal Specula

Integrated light source that is ready to use out of the package,  
simply pull switch when ready.



# Welch Allyn

## KleenSpec® Single Use LED Vaginal Specula



### LED Technology

- > Provides continuous output for 30 minutes
- > Emits crisp white light with no blue or yellow hues for true visualization of tissue color
- > Produces uniform light output, with no hot spots, to help provide the visualization needed for improved patient exams



- > LED light source is sealed within handle to reduce patient risk and placement provides unobstructed view of exam area
- > Smooth molded edges designed to provide maximum patient comfort
- > 100% Acrylic, as opposed to an acrylic-styrene mix, is 50% stronger minimizing the chance of cracking or breaking
- > Ready to use out of the package, simply pull switch when ready



- > Lithium-primary battery provides longer shelf life and withstands extreme temperature variations better than alkaline
- > Option to remove lithium-primary battery and LED for disposal or recycling
- > Designed for quality and durability from a trusted manufacturer and made in the U.S.A



590XS-LED	KleenSpec® Single Use LED Vaginal Specula, Extra Small (24/Box, 4 Boxes/Case)
59000-LED	KleenSpec® Single Use LED Vaginal Specula, Small (24/Box, 4 Boxes/Case)
59001-LED	KleenSpec® Single Use LED Vaginal Specula, Medium (24/Box, 4 Boxes/Case)

For more information, please contact your local Welch Allyn representative or call **1.800.535.6663**