

## Cordless Illumination System

Advanced LED light source provides improved visualization of the exam area by supplying brighter, whiter light.

### Cordless Illumination System

- No cord to get in the way
- Use with KleenSpec® Disposable Specula
- Automatic on/off upon specula insertion and removal



Learn more  
on page 188



## Specula Lighting Systems

KleenSpec® Cordless Illumination System .....	188
KleenSpec Corded Illumination System.....	188

## Vaginal Specula

KleenSpec Vaginal Specula .....	189
KleenSpec Vaginal Specula with Attached Sheath .....	190
KleenSpec Vaginal Specula with Smoke Tube .....	191
KleenSpec Economy Series .....	191

## Colposcope

Video Colposcope .....	192
Image Capture System.....	193

# SPECULA LIGHTING SYSTEMS


## WOMEN'S HEALTH



**79910**  
KleenSpec® Cordless  
Illuminator



**78810**  
KleenSpec® Corded  
Illumination System



Transformer  
with 5 ft Cord

**78800**  
Light Handle with On/Off Switch  
(shown with 73410, Transformer  
with 5 ft cord)

Halogen Lamp

**78820**  
KleenSpec® Corded  
Illumination Storage  
System



## KleenSpec® Cordless Illumination System

Designed with enhanced unique new features, including bright, white light with our new LED light sources. Single-use Sheath for KleenSpec Cordless Illumination Systems.

- Advanced LED light source provides improved visualization of the exam area by supplying brighter, whiter light
- No cord to get in the way during procedures
- Automatically illuminates upon insertion of the illuminator and once it is removed the light goes off
- Use your Cordless Illumination System as an all-purpose light source as well as with the KleenSpec® Disposable Vaginal Speculum

**79910** Complete KleenSpec Cordless Illuminator with Charging Station

**79920** Complete KleenSpec Cordless Illuminator with Charging Station and Auxiliary Power Cord

**79900** KleenSpec Cordless Illuminator

**73910** Charging Station for Cordless Illuminator

**75910** Auxiliary Power Cord Accessory for Cordless Illuminator

## KleenSpec® Cordless Illumination System: Accessories

**59010** Disposable Sheath for the Cordless Illumination System

## KleenSpec® Corded Illumination System

Designed with enhanced unique new features.

- Unique integrated light source transmits improved illumination for accurate diagnosis
- Integrated on/off switch close to the task area for operational convenience
- New ergonomic design allows for easy cleaning
- Reengineered, more robust cable for your busy practice needs
- New lamp design with improved lamp life

**78810** KleenSpec Corded Illumination System—Light Handle with Integrated On/Off Switch and Transformer

**78800** Light Handle with On/Off Switch Only

## KleenSpec Corded Illumination System: Accessories

**73410** Transformer with 5 ft Cord

**08800-U** 4.6 V Halogen Replacement Lamp for 78810 Illumination System

**78820** KleenSpec® Corded Illumination Storage System

### KleenSpec® Vaginal Specula

Welch Allyn's single-use specula eliminate the time and expense of reprocessing and reduce the chance of cross-contamination. Each speculum is designed and manufactured to deliver durability and reliability.

- Newly designed light pipe provides a homogeneous spot and eliminates glare-back
- New wider handle provides improved ergonomics, better balance, and is much easier to manipulate during exams
- Smooth, molded acrylic will not pinch or bind tissue for greater patient comfort during the exam
- The new 590 Series KleenSpec Vaginal Speculum is compatible with the new cordless illumination system (#79910) as well as the illumination system (#78810)
- Pederson Style Corded

**59000** KleenSpec 590 Series Disposable Vaginal Specula, Small (24/box; 4 boxes/case)

**59001** KleenSpec 590 Series Disposable Vaginal Specula, Medium(24/box; 4 boxes/case)

**59004** KleenSpec 590 Series Disposable Vaginal Specula, Large (19/box; 4 boxes/case)

### KleenSpec Vaginal Specula: Accessories

**58014** KleenSpec Vaginal Specula Storage Dispenser



**59000**  
590 Series Disposable  
Vaginal Specula, Small



**59001**  
590 Series Disposable  
Vaginal Specula, Medium



**59004**  
590 Series Disposable  
Vaginal Specula, Large



**58001S**  
580 Series Disposable  
Vaginal Specula with Sheath,  
Medium (shown with 78810,  
Illumination System)

## KleenSpec® Vaginal Specula with Attached Sheath

Reduces the risk of illuminator cross-contamination.

- Preattached sheath easily pulls down over the Welch Allyn corded illumination system
- The deployed sheath protects the light cord from cross-contamination
- Unique Welch Allyn high-quality molded acrylic for patient comfort and reliability
- Only compatible with corded illumination system (catalog #78810)
- Penderson Style

**58000S** KleenSpec 580 Series Disposable Vaginal Specula with Sheath, Small (20/box; 4 boxes/case)

**58001S** KleenSpec 580 Series Disposable Vaginal Specula with Sheath, Medium (20/box; 4 boxes/case)

**58004S** KleenSpec 580 Series Disposable Vaginal Specula with Sheath, Large (15/box; 4 boxes/case)



**58000S, 58001S and 58004S**  
580 Series Disposable Vaginal  
Specula with Sheath

### KleenSpec® Vaginal Specula with Smoke Tube

For smoke evacuation during electrosurgical procedures.

- Built-in smoke evacuation channel that assures a clear view
- Convenient single-use disposable for patient safety
- Excellent visibility—transparent blades allow full visualization of the exam area
- Smooth, molded acrylic will not pinch or bind tissue for greater patient comfort during exams
- The illumination system is completely sealed when inserted into the speculum
- Easy one-handed operation provides multiple settings with positive locking positions
- Pederson Style

**59005** KleenSpec 590 Series Disposable Vaginal Specula with Smoke Evacuator Tube, Small (11/box; 4 boxes/case)

**59006** KleenSpec 590 Series Disposable Vaginal Specula with Smoke Evacuator Tube, Medium (11/box; 4 boxes/case)



**59005**  
KleenSpec® 590 Series  
Disposable Vaginal  
Specula with Smoke  
Evacuator Tube, Small

### KleenSpec® Economy Series

Integrated illumination in an economical package.

- No sharp edges with this quality economy line
- Expendable illuminator can be re-sterilized up to 50 times before disposal
- No need for bulky exam lights
- Graves Style

**58600** KleenSpec 586 Series Disposable Vaginal Specula, Small (25/box; 4 boxes/case)

**58601** KleenSpec 586 Series Disposable Vaginal Specula, Medium (25/box; 4 boxes/case)



**58600**  
KleenSpec® 586 Series  
Disposable Vaginal  
Specula, Small with  
78800 light handle

### KleenSpec Economy Series: Accessories

**78600** Expendable Illuminator for Vaginal Specula (6/box, 6 boxes/case)



**88000A**  
Video Colposcope  
with Vertical Stand

## Video Colposcope

Offer more complete care with true video colposcopy combining light, power and video for high-resolution imaging, superior diagnostics and information management.

- The direct-to-screen colposcope gives you full-screen, high-resolution images with no couplers or beam splitters to degrade image quality or field of view
- Welch Allyn Solarc® lighting generates 50% whiter, brighter light for true tissue color
- Push-button on-screen magnification index
- Unique polarization filter removes glare for accurate tissue assessment
- Electronic green filter removes red from image with no light loss
- Store, retrieve, compare, annotate and transmit images with advanced image management
- A variety of documentation options are available

**88000A** Video Colposcope with Vertical Stand

**89000A** Video Colposcope with Swing Arm

**88001A** Video Colposcope with Vertical Stand and Monitor

**89001A** Video Colposcope with Swing Arm and Monitor

**88007** Video Colposcope with Vertical Stand, Monitor, Dust Cover and Replacement Lamp

## Video Colposcope: Accessories

**88040** Dust Cover for Vertical Stand Video Colposcope

**89040** Dust Cover for Swing Arm Video Colposcope

**09800-U** Replacement Solarc Lamp Assembly for Video Colposcope

**88010** Cervical Model

**65176** Sony® 14" LCD Monitor

## Printers & Accessories

**65189** Sony UP20 Color Video Printer

**88500** RS-232 Interface Cable (Replacement)

**88600** S-Video Cable (Replacement)


**65164** FS-20 Foot Switch (one function)

**65172** JVC® S-VHS Recorder Printers and Accessories

**65184** UPC 21S Print Pack—3 sets of 80 prints with ink cartridges (for UP20 printer model)



**89000A**  
Video Colposcope  
with Swing Arm



Welch Allyn Solarc lighting generates 50% whiter, brighter light for true tissue color

## Image Capture System

Take full advantage of your Welch Allyn Video Colposcope with this simple and affordable software-based imaging solution. It works directly with the colposcope to provide imaging documentation that can be integrated directly into your existing office technology and workflow.

- Capture, review, print and export colposcope images quickly and easily
- View captured images in real time—no more waiting for hard copies to print
- Electronic image storage facilitates patient tracking and management
- Communicates directly with the Welch Allyn Video Colposcope
- Flexible, affordable and easy to use
- Images are saved in standard file formats, making them easy to integrate with an existing patient management system

**88801** Image Capture System for Desktop PC (computer not included)

**88802** Image Capture System for Laptop PC (computer not included)



**88802**

Image Capture System for Laptop PC  
(computer not included)