



Nearly 100 years in the making...

Welch Allyn 767 Wall and Mobile Aneroids

Improved aesthetics and a new modern design make the 767 Aneroid an easy-to-use and attractive blood pressure aneroid. Based on Welch Allyn's near century of experience building and designing innovative blood pressure technology, it's the smart choice for accurate and reliable readings in your facility.





Benefits that matter.

The 767 Aneroid is a testament to Welch Allyn's nearly 100-year history of blood pressure innovation—combining the traditional benefits you have come to expect from Welch Allyn with improved aesthetics and a new, contemporary design.

Accurate

- Laser-engraved dial face provides unmatched accuracy
- Certified accurate to ±3 mm Hg
- Recessed dial face and highcontrast red pointer increase visibility and reduce parallax errors

Reliable

- Jeweled movement contributes to a longer life
- Lifetime calibration warranty proves our commitment to durability

Convenient

- New 5-leg stand provides stability and maneuverability
- Built-in storage basket on wall and desktop model

Safe

- FlexiPort® blood pressure cuffs are antimicrobial treated to help prevent mildew and bacteria growth
- Latex-free
- Mercury-free

Ordering Information

767 Wall Aneroid

7670-01 767 Wall Aneroid with Reusable One-Piece Adult Cuff*
7670-01C 767 Wall Aneroid with Reusable One-Piece Child Cuff

7670-01KPA 767 Wall Aneroid (KPA) with Adult Cuff

7670-02 767 Wall Aneroid with 8 ft. Coiled Tubing (No cuff included)

767 Mobile Aneroid

7670-03 767 Mobile Aneroid on 5-Leg Stand with Reusable One-Piece Adult Cuff*

7670-03K 767 Mobile Aneroid on 5-Leg Stand Only (No cuff included)
7670-10 767 Mobile Aneroid with Reusable One-Piece Adult Cuff*
7670-10K 767 Mobile Aneroid with Reusable One-Piece Child Cuff

7670-10KPA 767 Mobile Aneroid (KPA) on 5-Leg Stand with Reusable One-Piece Adult Cuff

7670-70 767 Aneroid with Mobile Bracket

767 Desktop Aneroid

7670-16 767 Desktop Model with Reusable One-Piece Adult Cuff

7670-16KPA 767 Desktop Model (KPA) with Adult Cuff

*Model available with traditional two-piece cuff and bladder inflations system. Add "CB" to the end of part number (e.g. 7670-01CB).

Not all models are available in all markets, check with your local Welch Allyn Representative. For more information about our blood pressure solutions, visit welchallyn.com.



Simple color changes and new rounded edges give the aneroid a modern look.





Advancing Frontline Care™

Welch Allyn Corporate Headquarters

4341 State Street Road, P.O. Box 220 Skaneateles Falls, NY 13153-0220 U.S.A. Telephone: +1.315.685.4602 +1.800.535.6663

Fax: +1.315.685.0054 E-mail: info@welchallyn.com Welch Allyn European Customer Service: +353 46 90 67790

Australia: +61 29 638 3000 China: +86 21 6327 9631 France: +33 155 69 58 49 Germany: +49 6950 985 132

Latin America (Miami office): +1 305 669 9003

Malaysia: +603 78753341 Netherlands: +31 202 061 360 Singapore: +65 6419 8100 South Africa: +27 11 777 7555 Sweden: +46 858 536 551 United Kingdom: +44 207 365 6780

©2010 Welch Allyn SM4022 Rev A