

seca 354

Digital baby scales,
also ideal as floor scales
for children



- Baby scale and floor scale for children in one.
- Weighing tray easy to detach from base.
- Breast-Milk-Intake-Function.
- Mobile and robust for transport.

seca 354:

Convertible and robust –
compact quality with dual purpose.

At home, in surgery and around the world wherever life begins: reliable weight control is decisive for a healthy start in life. The seca 354 is an all-around scale which has been developed on the basis of practical experience. It can be used not only as a baby scale but also as a floor scale for small children.



Baby scale and floor scale for children in one.

Babies can be weighed safely and comfortably in the tray of the seca 354. The rounded shape of the tray with its slightly raised sides gives the baby a feeling of security. By detaching the tray from the base, the scale is converted into a fully functional digital floor scale with a capacity of 44 lbs.

Very safe due to secure attachment of tray.

In no time at all, the baby scale can be converted into a floor scale: pressing a button unfastens the mechanism attaching the tray to the base. When locked, this secure attachment ensures that the baby scale is easy to lift and safe to transport.

Easy to operate due to self-explanatory button functions.

The seca 354 is easy to operate due to its intuitive operating concept with self-explanatory buttons. The auto-HOLD function means that, after being weighed, the baby can first be taken care of before one has to read off the result. With the TARE function, the weight of nappies or soft mats can be reliably deducted. With the Breast-Milk-Intake-Function, one can determine how much milk the baby has drunk when breast fed. Because, at the touch of a button, the weight of the baby determined before feeding is deducted from the weight determined after feeding.



Mobile and robust for transport.

Its compact, robust construction predestines the seca 354 as the scale for all those who have to be mobile in the interest of the health of their youngest patients. The power-saving battery operation – up to 20,000 weighings with just one set of batteries – makes the scale independent of electricity power supply.



seca 414



seca 413



The weighing tray can be removed quickly and easily at the touch of a button.



seca 354

Technical Data

- Capacity: 44 lbs / 20 kg
- Graduation: 0.5 oz < 22 lbs > 1 oz, 10 g < 10 kg > 20 g
- Dimensions (WxHxD): 21.7 x 6.1 x 13.1" / 552 x 156 x 332 mm
- Dimensions/tray (WxHxD): 20.7 x 3.1 x 9.8" / 525 x 80 x 250 mm
- Dimensions/platform (WxHxD): 11 x 2.6 x 9.3" / 280 x 65 x 235 mm
- Weight: 5 lbs / 2.3 kg
- Power supply: batteries
- Functions: TARE, BMIF, auto-HOLD, lbs/kg switch-over, automatic switch-off
- Optional: carrying case seca 414 and seca 413