


# NEW! HEINE® iC 1

## 21<sup>st</sup> Century Digital Dermatoscopy



# NEW! HEINE® iC1 – 21<sup>st</sup> Century Digital Dermatology

The HEINE iC1 is a new image capturing device to record clinical and dermatoscopic images for documenting, monitoring and supporting differential diagnosis of melanoma and other skin lesions. Innovative HEINE LED<sup>HQ</sup> illumination combined with precision optics provide high resolution and reflex-free images with an Apple® iPhone® 5, 5s, 6, 6s and SE. The integrated Toggle Function allows you to easily alternate between polarised and non-polarised views.



## [ ACCESSORIES ]

<b>USB cord with E4-USB Medical approved plug-in power supply</b> [01] + [02] Medical approved plug-in power supply with 5V and 500mA; changeable plug adapter	Art.No: <b>X-000.99.303</b>
<b>Contact plate</b> with scale [03]	Art.No: <b>K-000.34.201</b>
<b>Mobile phone case</b> , for Apple iPhone 6/6s [04]	Art.No: <b>K-000.34.250</b>
<b>Mobile phone case</b> , for Apple iPhone 5/5s and iPhone SE	Art.No: <b>K-000.34.251</b>



## [ TECHNICAL DATA ]

Catalogue number: <b>K-270.28.305</b> <b>iC1 set/6</b>	complete with: HEINE iC1 with contact plate with scale, mobile phone case iC1/6 currently available for Apple iPhone 6/6s [05], USB cord with medical approved plug-in power supply
Catalogue number: <b>K-271.28.305</b> <b>iC1 set/5</b>	complete with: HEINE iC1 with contact plate with scale, mobile phone case iC1/5 currently available for Apple iPhone 5/5s and iPhone SE, USB cord with medical approved plug-in power supply
Available versions	iC1/5 for iPhone 5/5s and iPhone SE; iC1/6 for iPhone 6/6s
Weight	140g
Dimensions	140x65x42 mm up to 155x75x50 mm
Users manual	DE, GB, FR, ES, IT, SE, NL, DK, FI, PT
Operating conditions	+10 °C to +35 °C, 30% to 75% RH, 700hPa to 1060hPa
Terms of transport	-20 °C to +50 °C, 45% to 85% RH, 500hPa to 1060hPa
Electrical data	micro USB (5V) <500mA, charging time typ. 2.5h, operating time typ. 3.5h Li-Polymer rechargeable battery
Optical data	Optical magnification = 15x, digital magnification = 40x (with iPhone 6 and 24" display), 21 mm field of view
LED	4x high power LEDs, colour temp. typ. 5000K, working life typ. 100000h, typ. CRI ≥80

\*The Apple iPhone 5/5s/6 or 6s is not included in the scope of delivery.  
Apple, the Apple logo and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries.  
App Store is a service mark of Apple Inc.

